

**ThermoFisher**  
SCIENTIFIC

# PCR solutions

Products designed for what you aspire to achieve



# Amplification made for breakthroughs

In your pursuit to advance science, every experiment matters. There's no time to start again, no need to wonder if the products you chose will set you back or propel you forward. With our comprehensive portfolio of thermal cyclers, PCR plastics, reagents, and service plans, you know exactly what you're getting—quality made to support your every win, big and small, at every step.

 Find out more at [thermofisher.com/pcrworkflow](https://www.thermofisher.com/pcrworkflow)

# Contents

Which instrument fits your needs?	4
ProFlex PCR System	6
VeritiPro Thermal Cycler	8
SimpliAmp Thermal Cycler	10
MiniAmp Thermal Cycler	12
Automated Thermal Cycler	14
Thermal Cycler Fleet Control Software	16
Thermal cycler services and support	16
MicroAmp PCR plastics	18
SuperScript IV Reverse Transcriptase	22
PCR enzymes	24
PCR primers and dNTPs	26
Nucleic acid separation and analysis	27
Specialty enzymes	31
Molecular biology education	31

# Which instrument fits your needs?

Providing precise, consistent results, Applied Biosystems™ thermal cyclers are designed for every challenge, application, and budget. Be confident in your results with instruments designed for reliability, accuracy, and intuitive programming. Reduce time spent on PCR optimizations with Applied Biosystems™ VeriFlex™ Block temperature control technology. Manage your instruments remotely with secure access.



	Ultimate flexibility and throughput	Ultimate performance
	ProFlex PCR System	VeritiPro Thermal Cycler
Max. sample throughput	480,000 reactions	384 reactions
Max. block ramp rate	6.0°C/sec	6.0°C/sec
Block formats (temperature optimization)	<ul style="list-style-type: none"> <li>• 3 x 32-well, 0.2 mL (2-zone VeriFlex Block)</li> <li>• 96-well, 0.2 mL (6-zone VeriFlex Block)</li> <li>• 2 x 96-well, 0.2 mL</li> <li>• 2 x flat block</li> <li>• 2 x 384-well, 0.02 mL</li> </ul>	<ul style="list-style-type: none"> <li>• 96-well, 0.2 mL (6-zone VeriFlex Block)</li> <li>• 384-well, 0.02 mL</li> </ul>

For Research Use Only. Not for use in diagnostic procedures.



= cloud-enabled instrument



= automation-ready instrument






Elegantly simple and precise	Routine PCR	Designed for easy robotic integration
SimpliAmp Thermal Cycler	MiniAmp Thermal Cycler	Automated Thermal Cycler
96 reactions	96 reactions	384 reactions
4.0°C/sec	3.0°C/sec	3.5°C/sec
<ul style="list-style-type: none"> <li>96-well, 0.2 mL (3-zone VeriFlex Block)</li> </ul>	<ul style="list-style-type: none"> <li>96-well, 0.2 mL</li> </ul>	<ul style="list-style-type: none"> <li>96-well, 0.2 mL, compatible with full- or semi-skirted plates</li> <li>384-well, 0.02 mL</li> </ul>

Looking for a thermal cycler for *in vitro* diagnostic use? Learn more about the Applied Biosystems™ Veriti™ Dx Thermal Cycler at [thermofisher.com/veritidx](http://thermofisher.com/veritidx)



**Don't forget reagents**—choose from industry-leading enzymes such as Invitrogen™ SuperScript™ IV reverse transcriptases and Invitrogen™ Platinum™ SuperFi™ II reagents (see pages 22–25).



Interested in private-label thermal cyclers or PCR plastics? To find out more, go to [thermofisher.com/oem](http://thermofisher.com/oem)

# Ultimate flexibility and throughput

## ProFlex PCR System

The Applied Biosystems™ ProFlex™ PCR System combines flexible configuration and control features to fit how you work today and tomorrow with the reliability you've come to expect from Applied Biosystems™ products. Interchangeable block formats allow you to maximize your throughput or run independent experiments concurrently.

The ProFlex PCR System is cloud-enabled, giving you the freedom to design and securely upload your methods, monitor runs, and check instrument availability from any mobile device or desktop computer with the Thermo Fisher™ Connect Platform.

- **Multi-user accessible**—run three experiments at once
- **Flexible block configuration**—accepts five different block formats for optimization and throughput
- **Cloud-enabled**—conveniently access your instrument anytime, anywhere with the Connect Platform
- **Fleet control-compatible**—manage multiple instruments, users, and methods with Applied Biosystems™ Thermal Cycler Fleet Control Software



 = cloud-enabled instrument

### Five interchangeable block options

The ProFlex PCR System has five different blocks that can be changed with the flip of a switch, including a 3 x 32-well block. This allows up to three experiments to be run simultaneously and completely independently of each other.

Dual 96-well and dual 384-well blocks are available for high-throughput needs. A dual flat block is also available to support Applied Biosystems™ OpenArray™ plate technology for genotyping analysis on the Applied Biosystems™ QuantStudio™ 12K Flex Real-Time PCR System as well as our sealed-chip technology on the QuantStudio™ 3D Digital PCR System.



3 x 32-well



96-well



Dual 96-well



Dual 384-well



Dual flat

 Find out more at [thermofisher.com/proflex](https://thermofisher.com/proflex)

## Specifications

Block format	3 x 32-well, 0.2 mL, 2-zone VeriFlex Block independent control	96-well, 0.2 mL, 6-zone VeriFlex Block	2 x 96-well, 0.2 mL	2 x flat block for chips and arrays	2 x 384-well, 0.02 mL
Features	Run three experiments at once or at different times	Perform complete optimization work with full 96-well VeriFlex Block	High throughput in 96-well format	High throughput capability: 8 x 3,072 OpenArray Plate* or 24 x 20K Chip**	High throughput in 384-well format
Max. block ramp rate	6.0°C/sec		3.0°C/sec	1.6°C/sec	3.0°C/sec
Max. sample ramp rate	4.4°C/sec		1.6°C/sec	NA	1.6°C/sec
Temperature accuracy	±0.25°C (35–99.9°C)				
Temperature range	0–100.0°C				
Temperature uniformity	<0.5°C (30 sec after reaching 95°C)				
Dimensions (H x W x D)	27.2 x 33.0 x 56.5 cm (10.6 x 13 x 22 in.)				
Weight	18,8 kg (41 lb)		20,4 kg (45 lb)		
PCR volume range	10–80 µL		10–100 µL	33 nL	5–20 µL
Instrument memory	USB, onboard				
Display interface	8.4-inch color TFT LCD				
Power	100–240 V, 50–60 Hz, max. 950 VA				
VeriFlex Block	2 temperature zones per block (5°C zone-to-zone)	6 temperature zones, 25°C range (5°C zone-to-zone)	NA		
Data connectivity†	Cloud or mobile via Ethernet or Wi-Fi				

\* OpenArray Plate is compatible with the QuantStudio 12K Flex Real-Time PCR System.

\*\* The Applied Biosystems™ QuantStudio™ 3D Digital PCR 20K Chip is compatible with the QuantStudio 3D Digital PCR System.

† The Instrument Connect app, available at Apple™ and Google™ app stores, can be used to monitor your instrument.

The Connect Platform can be used to create and share protocols, schedule an instrument, and start and monitor runs remotely.

## Ordering information

Product	Complete system Cat. No.	Block only Cat. No.	Instrument + 5-year* Rapid Exchange warranty Cat. No.
ProFlex 96-Well PCR System	4484075	4483637	A27934
ProFlex 3 x 32-Well PCR System	4484073	4483638	A28986
ProFlex Dual 96-Well PCR System	4484076	4484071	A27937
ProFlex Dual Flat PCR System	4484078	4484074	A27931
ProFlex Dual 384-Well PCR System	4484077	4484072	A30229

## Recommended plastics

3 x 32-well block	Cat. No.	96-well block	Cat. No.	384-well block	Cat. No.
MicroAmp TriFlex 3 x 32-Well PCR Reaction Plate	A32811	MicroAmp EnduraPlate Optical 96-Well Clear Reaction Plate with Barcode	4483354	MicroAmp Optical 384-Well Reaction Plate	4343370
MicroAmp 8-Tube Strip with Attached Domed Caps, 0.2 mL	A30589	MicroAmp TriFlex 3 x 32-Well PCR Reaction Plate	A32811	MicroAmp Optical 384-Well Reaction Plate with Barcode	4309849
MicroAmp Reaction Tube with Cap, 0.2 mL	N8010540	MicroAmp 8-Tube Strip with Attached Domed Caps, 0.2 mL	A30589	MicroAmp EnduraPlate Optical 384-Well Clear Reaction Plate with Barcode	4483285
MicroAmp 32-Well Clear Adhesive Film	A32812	MicroAmp Clear Adhesive Film	4306311	MicroAmp Clear Adhesive Film	4306311
		MicroAmp 32-Well Clear Adhesive Film	A32812		

\* The system comes with a 2-year standard warranty; 3 years of an extended warranty may be added for a total of 5 years of service.



## Did you know?

The ProFlex PCR System, Applied Biosystems™ VeritiPro™ Thermal Cycler, and the Applied Biosystems™ SimpliAmp™ Thermal Cycler feature VeriFlex temperature control technology, which enables more precise and efficient PCR optimization.

Find out more at [thermofisher.com/veriflextechnology](https://thermofisher.com/veriflextechnology)


# Ultimate performance

## VeritiPro Thermal Cycler

The VeritiPro Thermal Cycler delivers proven reliability with advanced temperature control technology and connectivity. Take advantage of next-level PCR optimization with the precision offered by VeriFlex temperature control technology. Connect to the cloud-enabled VeritiPro Thermal Cycler remotely. Conveniently design and securely upload your methods, monitor runs, and check instrument availability from any mobile device or desktop computer with the Connect Platform.

- **Innovative design**—ramp rate of 6.0°C/sec, quiet fan, and ergonomic soft-close lid
- **VeriFlex Block**—precise control over six independent temperature zones for PCR optimization
- **Fleet control-compatible**—manage multiple instruments, users, and methods with Thermal Cycler Fleet Control Software




 = cloud-enabled instrument

### Precise control with VeriFlex temperature control technology

VeriFlex Blocks are constructed of segmented metal blocks with separate heating and cooling elements below each, enabling:

- **More precise control over PCR optimization**—each block can be set with up to six specific temperatures
- **Precise incubation**—use the six temperature zones of the VeriFlex Block to do enzyme studies, restriction digests, or any other process that requires precise temperature control



 Find out more at [thermofisher.com/veritipro](https://thermofisher.com/veritipro)



### Tired of water baths?

Incubate samples at up to six different temperatures simultaneously for enzyme activation studies, restriction digests, or sequencing library preps with the VeritiPro Thermal Cycler.



## Specifications

Block format	96-well, 0.2 mL, 6-zone VeriFlex Block	384-well, 0.02 mL, isothermal
Max. block ramp rate	6.0°C/sec	5.0°C/sec
Max. sample ramp rate	4.4°C/sec	3.5°C/sec
Temperature accuracy	±0.25°C (35–99.9°C)	
Temperature range	0–100.0°C	
Temperature uniformity	<0.5°C (30 sec after reaching 95°C)	
Dimensions (H x W x D)	21.7 x 24.5 x 46.5 cm (8.5 x 9.6 x 18.3 in.)	
Weight	12.0 kg (26.5 lb)	
PCR volume range	10–100 µL	5–20 µL
Instrument memory	USB port and 16 GB onboard memory; onboard capacity >1,000 protocols	
Display interface	8-inch color TFT LCD	
Power	100–240 V, 50–60 Hz, max. 700 W	
VeriFlex Block range	30°C range across block, 6 temperature zones (up to 10°C zone-to-zone)	NA
Data connectivity*	Cloud or mobile via Ethernet or Wi-Fi	

\* The Instrument Connect app, available at Apple and Google app stores, can be used to monitor your instrument.

The Connect Platform can be used to create and share protocols, schedule an instrument, and start and monitor runs remotely.

## Ordering information

Product	Cat. No.	Instrument + 5-year* Rapid Exchange warranty Cat. No.
VeritiPro Thermal Cycler, 384 well	A48140	A50326
VeritiPro Thermal Cycler, 96 well	A48141	A48765

### Recommended plastics

384-well block	Cat. No.
MicroAmp Optical 384-Well Reaction Plate	4343370
MicroAmp Optical 384-Well Reaction Plate with Barcode	4309849
MicroAmp EnduraPlate Optical 384-Well Clear Reaction Plate with Barcode	4483285
MicroAmp Clear Adhesive Film	4306311
96-well block	Cat. No.
MicroAmp EnduraPlate Optical 96-Well Clear Reaction Plate with Barcode	4483354
MicroAmp Optical 96-Well Reaction Plate	N8010560
MicroAmp Clear Adhesive Film	4306311
MicroAmp TriFlex 3 x 32-Well PCR Reaction Plate	A32811
MicroAmp 32-Well Clear Adhesive Film	A32812
MicroAmp 8-Tube Strip with Attached Domed Caps, 0.2 mL	A30589

\* The system comes with a 2-year standard warranty; 3 years of an extended warranty may be added for a total of 5 years of service.

# Elegantly simple and precise

## SimpliAmp Thermal Cycler

The SimpliAmp Thermal Cycler is an easy-to-use, compact, and accurate thermal cycler designed to fit every lab's essential PCR workflow. Features like a responsive color touchscreen and VeriFlex temperature control technology enable simple, accurate optimization. Plus, the SimpliAmp Thermal Cycler is cloud-enabled, giving you the freedom to design and securely upload your methods, monitor runs, and check instrument availability from any mobile device or desktop computer with the Connect Platform.

- **Intuitive interface**—large, easy-to-use color touchscreen for easy programming and quick status checks
- **VeriFlex Block**—three independent temperature zones for PCR optimization
- **Cloud-enabled**—conveniently access your instrument anytime, anywhere with the Connect Platform
- **Compact design**—helps save bench space
- **Fleet control-compatible**—manage multiple instruments, users, and methods with Thermal Cycler Fleet Control Software




 = cloud-enabled instrument

### Simple, intuitive user interface

The SimpliAmp Thermal Cycler has an 8-inch color touchscreen, making navigation of the intuitive menu options fast and efficient.

The status dial displays the current block temperature and elapsed run time.



Create a new run method using default templates or existing methods.

Select one of your existing run methods to start a run.

Use the instrument as a precise incubator for non-PCR workflows.

Find out more at [thermofisher.com/simpliamp](https://thermofisher.com/simpliamp)

## Specifications

<b>Block format</b>	96-well, 0.2 mL, 3-zone VeriFlex Block	
<b>Features</b>	<ul style="list-style-type: none"> <li>Enabled to run Fast chemistry</li> <li>Controllable ramp rate</li> <li>Program overwrite protection</li> </ul>	<ul style="list-style-type: none"> <li>Auto restart (after power outages)</li> <li>Edit program during experiment</li> <li>One-touch incubation</li> </ul>
<b>Max. block ramp rate*</b>	4°C/sec	
<b>Max. sample ramp rate*</b>	3°C/sec	
<b>Temperature accuracy</b>	±0.25°C (35–99.9°C)	
<b>Temperature range</b>	0–100.0°C	
<b>Temperature uniformity</b>	<0.5°C (30 sec after reaching 95°C)	
<b>Temperature calibration</b>	Calibrated to standards traceable to the National Institute of Standards and Technology (NIST)	
<b>Dimensions (H x W x D)</b>	21.0 x 24.0 x 46.0 cm (8.3 x 9.5 x 18.1 in.)	
<b>Weight</b>	8.3 kg (18.3 lb)	
<b>PCR volume range</b>	10–100 µL	
<b>Instrument memory</b>	2,000 MB onboard memory (capacity for >1,000 protocols); USB port for additional external storage	
<b>Display interface</b>	8-inch color TFT LCD	
<b>Power</b>	100–240 V, 50–60 Hz, max. 600 W	
<b>VeriFlex Block</b>	3 temperature zones, 20°C range (10°C zone-to-zone)	
<b>Data connectivity**</b>	Cloud or mobile via Ethernet or Wi-Fi	

\* At a reaction volume of 1 µL.

\*\* The Instrument Connect app, available at Apple and Google app stores, can be used to monitor your instrument. The Connect Platform can be used to create and share protocols, schedule an instrument, and start and monitor runs remotely.

## Ordering information

<b>Product</b>	<b>Cat. No.</b>	<b>Instrument + 5-year* warranty with Rapid Exchange plan Cat. No.</b>
SimpliAmp Thermal Cycler	A24811	A27603
High-Power USB Wi-Fi Module	A26774	NA

<b>Recommended plastics</b>	
<b>96-well block</b>	<b>Cat. No.</b>
MicroAmp EnduraPlate Optical 96-Well Clear Reaction Plate with Barcode	4483354
MicroAmp Optical 96-Well Reaction Plate	N8010560
MicroAmp Clear Adhesive Film	4306311
MicroAmp TriFlex 3 x 32-Well PCR Reaction Plate	A32811
MicroAmp 32-Well Clear Adhesive Film	A32812
MicroAmp 8-Tube Strip with Attached Domed Caps, 0.2 mL	A30589



## Did you know?

ProFlex, VeritiPro, SimpliAmp, and Applied Biosystems™ MiniAmp™ Thermal Cyclers have thermal simulation modes that make the transition from other thermal cyclers simple, accurate, and efficient. A library of modes that mimic the ramp rates of other instruments is available on each instrument.


# Routine PCR

## MiniAmp Thermal Cycler

The MiniAmp Thermal Cycler delivers the reliability you've come to expect from Applied Biosystems™ technology for routine PCR at every lab bench.

- **Compact design**—fits in everyone's lab space at just 7.5 in. (19 cm) wide
- **Cloud-enabled**—conveniently access your instrument anytime, anywhere with the Connect Platform
- **Fleet control-compatible**—manage multiple instruments, users, and methods with Thermal Cycler Fleet Control Software




 = cloud-enabled instrument

### Secure remote access with cloud-enabled instruments

The ProFlex, VeritiPro, SimpliAmp, and MiniAmp Thermal Cyclers are cloud-enabled instruments, allowing you secure, private access with an account on the Connect Platform. Anywhere, anytime with any mobile device or desktop computer you can:

- Design and share protocols
- Schedule an instrument
- Start or stop a run
- Check run status



 Find out more at [thermofisher.com/miniamp](https://thermofisher.com/miniamp)



### Did you know?

Most Applied Biosystems thermal cyclers come with a two-year standard warranty and a starter kit that includes Applied Biosystems™ PCR plate and tube samples and all the tools you need. Packages that include the instrument, an extended warranty, and consumables are available. Visit the instrument webpage to view available packages.



## Specifications

	MiniAmp Thermal Cycler	MiniAmp Plus Thermal Cycler
Block format	96-well, 0.2 mL isothermal block	96-well, 0.2 mL, 3-zone VeriFlex Block
Max. block ramp rate*	3.0°C/sec	3.5°C/sec
Max. sample ramp rate*	2.2°C/sec	2.7°C/sec
Temperature accuracy	±0.25°C (35–99.9°C)	
Temperature range	0–100.0°C	
Temperature uniformity	<0.5°C (30 sec after reaching 95°C)	
Temperature calibration	Calibrated to standards traceable to the NIST	
Dimensions (H x W x D)	20 x 19 x 39 cm (7.9 x 7.5 x 15.4 in.)	
Weight	5.9 kg (13.0 lb)	
PCR volume range	10–100 µL	
Instrument memory	2,000 MB onboard memory (capacity for >1,000 protocols); USB port for additional external storage	
Display interface	5-inch color TFT LCD	
Power	100–240 V, 50–60 Hz, max. 500 W	
VeriFlex Block	NA	3 temperature zones, 20°C range (10°C zone-to-zone)
Data connectivity**	Cloud or mobile via Ethernet or Wi-Fi	

\* At a reaction volume of 1 µL.

\*\* The Instrument Connect app, available at Apple and Google app stores, can be used to monitor your instrument.

The Connect Platform can be used to create and share protocols, schedule an instrument, and start and monitor runs remotely.

## Ordering information

Product	Cat. No.	Instrument + 5-year warranty with Rapid Exchange plan Cat. No.
MiniAmp Thermal Cycler	A37834	A38081
MiniAmp Plus Thermal Cycler	A37835	A38077
High-Power USB Wi-Fi Module	A26774	NA

Recommended plastics	
96-well block	Cat. No.
MicroAmp EnduraPlate Optical 96-Well Clear Reaction Plate with Barcode	4483354
MicroAmp Optical 96-Well Reaction Plate	N8010560
MicroAmp Clear Adhesive Film	4306311
MicroAmp TriFlex 3 x 32-Well PCR Reaction Plate	A32811
MicroAmp 32-Well Clear Adhesive Film	A32812
MicroAmp 8-Tube Strip with Attached Domed Caps, 0.2 mL	A30589




Want to learn more about thermal cycler ramp rates and how they are calculated?  
Go to [thermofisher.com/ramprate](https://thermofisher.com/ramprate)

# Designed for easy robotic integration

## Automated Thermal Cycler

The Applied Biosystems™ Automated Thermal Cycler offers the flexibility, reliability, and performance needed in a complete PCR automation system. The small, easy-to-integrate format of the Automated Thermal Cycler (ATC) enables hands-free PCR results.

- **Flexible modular design and small footprint**—helps save space on deck
- **Automated lid**—easy, hands-free operation with a liquid handler or plate stacker
- **Free desktop software**—for PCR optimization prior to robotic integration
- **Plug-and-play drivers and Standardization in Lab Automation (SiLA) compatibility**—for easy integration on your liquid handler of choice
- **Fleet control-compatible**—manage multiple instruments, users, and methods with Thermal Cycler Fleet Control Software



= automation-ready instrument

## Designed for any stage of your workflow automation journey



### Stand-alone

Optimize assays before robotic integration with our direct software.



### Plug-and-play drivers

Ask about available drivers for leading robotic platforms.



### SiLA-compatible

Maximum robotic platform integration flexibility with SiLA rapid integration coding.

Find out more at [thermofisher.com/atc](https://thermofisher.com/atc)

## Specifications

Block format	96-well, 0.2 mL*	384-well, 0.02 mL
PCR volume range	10–100 $\mu$ L for full-skirted plates; 20–100 $\mu$ L for semi-skirted plates	5–20 $\mu$ L for full-skirted plates
Hardware integration features	<ul style="list-style-type: none"> <li>• Predrilled mounting and/or alignment points at each corner of the chassis</li> <li>• 3-side and top-plate access</li> <li>• Available in 3-connector configurations</li> </ul>	
Software features	<ul style="list-style-type: none"> <li>• Application programming interfaces (APIs) available for integration with robotics systems</li> <li>• SiLA Rapid Integration software—standardized programming access**</li> <li>• Free software available for instrument demonstration and stand-alone operation</li> </ul>	
Block module dimensions (H x W x D)	13.3 x 17.9 x 31.7 cm (5.2 x 7.0 x 12.5 in.)	
Control module dimensions (H x W x D)	7.0 x 25.7 x 33.1 cm (2.8 x 10.1 x 13.0 in.)	
Temperature accuracy	$\pm$ 0.25°C (35.0–99.9°C)	
Max. block ramp rate†	3.5°C/sec	2.8°C/sec
Max. sample ramp rate†	1.8°C/sec	1.6°C/sec
Temperature range	4–105°C (no condensation risk with subambient temperatures)	
Temperature uniformity	$\leq$ 0.50°C (20 sec after reaching 95°C)	
Temperature calibration	Calibrated to standards traceable to the NIST	
Service options	<ul style="list-style-type: none"> <li>• 2-year standard warranty includes Rapid Exchange service plan</li> <li>• NIST-traceable temperature probe equipment available</li> </ul>	
Weight	9.4 kg (20.7 lb) total (block module 6.0 kg (13.2 lb), control module 3.4 kg (7.5 lb))	
Power	100–240 V, 50–60 Hz, max. 600 W	
Flexible ramp rates	Program your own ramp rates, or use preprogrammed simulation modes	
Data connectivity	LAN	

\* Compatible with full- or semi-skirted plates. ATC semi-skirted adaptor required for use with semi-skirted 96-well plates. ATC ships with the full-skirted adaptor installed, which is required for use with full-skirted 96-well plates.

\*\* [sila-standard.org](http://sila-standard.org)

† At a reaction volume of 1  $\mu$ L.

## Ordering information

Product	96-well Cat. No.	384-well Cat. No.	3-year extended warranty* with Rapid Exchange plan Cat. No.
Automated Thermal Cycler System, laptop, 1 m cable	A31486	A33977	ZGEXSCATC3Y/ZGEXSCATC3843Y
Automated Thermal Cycler System, laptop, 3 m cable	A31487	A33978	
Automated Thermal Cycler System, laptop, 10 cm cable	A31488	A33979	
Automated Thermal Cycler System, 1 m cable	A31489	A33980	
Automated Thermal Cycler System, 3 m cable	A31490	A33981	
Automated Thermal Cycler System, 10 cm cable	A31491	A33982	
Automated Thermal Cycler Semi-Skirted Adaptor	A33044	NA	NA
Automated Thermal Cycler Full-Skirted Adaptor	A33045	NA	NA

\* The Automated Thermal Cycler comes standard with a 27-month warranty.

Recommended plastics			
96-well block	Cat. No.	384-well block	Cat. No.
MicroAmp EnduraPlate Optical 96-Well Full-Skirted Plate with Barcode, Clear*	A31728	MicroAmp EnduraPlate Optical 384-Well Clear Reaction Plate with Barcode	4483273
MicroAmp EnduraPlate Optical 96-Well Multicolor Reaction Plate with Barcode, semi-skirted**	4483356	MicroAmp EnduraPlate Optical 384-Well Multicolor Reaction Plate with Barcode	4483317
MicroAmp Clear Adhesive Film	4306311	MicroAmp Clear Adhesive Film	4306311

\* Available in multiple colors.

\*\* Requires semi-skirted plate adaptor (Cat. No. A33044); included with instrument purchase.



The Applied Biosystems™ MicroAmp™ EnduraPlate™ plate (Cat. No. A31728) was specifically designed for the Automated Thermal Cycler and has three-way barcoding. Learn about the amplification uniformity benefits at [thermofisher.com/automatedpcrplate](http://thermofisher.com/automatedpcrplate)

# Manage multiple users, methods, and instruments with a single intuitive interface

## Thermal Cycler Fleet Control Software



The Thermal Cycler Fleet Control Software is as intuitive and easy to use as our thermal cyclers are. With this tool, you'll be able to view and control all Applied Biosystems™ thermal cyclers in your facility.

**Secure**—comes preinstalled on a server that sits securely behind your firewall

**Flexible**—create custom permissions and rules for users, methods, and instruments

**Powerful**—quickly view your thermal cycler inventory and availability right from the dashboard

**Efficient**—track system changes with an audit log, and identify methods or instruments via barcode data entry

**Convenient**—control your thermal cyclers, reports, methods, and firmware upgrades from your office, another lab computer, or any other networked device

Find out more at [thermofisher.com/fleetcontrol](http://thermofisher.com/fleetcontrol)


# Science can't wait on instrument downtime

## Thermal cycler services and support

You can't afford downtime in your lab. Built on more than 35 years of service expertise, our superior service solutions for Applied Biosystems™ instruments and applications help keep your lab up and running. More than 1,400 trained professionals make up the industry's largest network, ready to assist you when you need it.

### Rapid Exchange

When fast turnaround time is critical and asset tag consistency is not required, choose our Rapid Exchange plan. If our Remote Service Center engineers can't repair your instrument through telephone support, we will ship you a factory-certified refurbished replacement instrument within one business day—and it's yours to keep. AB Repair Center (ABRC) Support Plus



### Off-site service plans at a glance

	Rapid Exchange	AB Repair Center (ABRC) Support Plus
Plan designed for	Labs that operate under critical timelines	Regulated labs that must maintain asset tag consistency to comply with regulated protocols
Repair time	Replacement shipped in 1 business day*	Loaner instrument shipped in 1 business day**
Off-site repair service including shipping, parts, and labor	✓	✓
Telephone and email access to our Remote Service Center	✓	✓

\* Subject to regional availability.

\*\* Your instrument repaired and returned in ~3 weeks. Service plan without loaner instrument also available.



Our AB Repair Center Support Plus plan is the best choice for customers in regulated environments that must maintain asset tag consistency to comply with regulated protocols and need their instrument returned after repair. With this mail-in option, we provide a loaner instrument for you to use while factory-trained engineers maintain and repair your instrument at our Remote Service Center. All instruments undergo certification performance tests, temperature verification, instrument and optical calibration, ground continuity, and firmware upgrades, when applicable.

Just mail us your instrument and we'll repair and ship it back to you within three weeks. All parts, labor, and shipping charges are included. Planned maintenance, temperature verification, and loaner instrument service\* can be added to any off-site service plan.

Check your instrument's coverage status at [thermofisher.com/manage](https://www.thermofisher.com/manage)

### Digital service innovations

Our instrument service plans include digital innovations that help keep your instruments and your lab running smoothly. With pioneering on-demand tools and capabilities such as remote support using mobile augmented-reality (AR) technology, instrument-driven support, and on-demand instrument training, we're constantly looking ahead so your lab never falls behind.

Explore our digital service innovations at [thermofisher.com/innovations](https://www.thermofisher.com/innovations)

### Temperature verification services

Multiplex dynamic temperature verification (MDTV) services are designed to deliver fast, accurate temperature readings and documentation to assess the accuracy and uniformity of your PCR systems. MDTV is available for Applied Biosystems thermal cyclers.\*\*

Have confidence in the accuracy of your results, learn more about the MDTV service at [thermofisher.com/mdtv](https://www.thermofisher.com/mdtv)

\* Loaner instruments are not available in all regions. Contact your services and support representative for availability.

\*\* MDTV service not currently offered on the 60-well Applied Biosystems™ GeneAmp™ PCR System 9700.

### Education services

It can be difficult to prepare yourself for what's next while you're focused on the work you have now. Our interactive professional training courses make it easier. We offer a combination of virtual and in-person classroom instruction, and hands-on learning in your lab to match your schedule, budget, and learning preferences. Whichever course style you choose, you'll learn from one of our 300 highly skilled application scientists who are available to lead sessions online, at your location, or at one of our 12 training centers located worldwide.

Explore courses at [thermofisher.com/educationservices](https://www.thermofisher.com/educationservices)


### Technical support

If you have questions about product selection or use, assay or experimental design, data analysis, or troubleshooting, contact our team of technical support scientists or access our online product and application support tools.

### How to reach us

To find your local support or technical support team, go to [thermofisher.com/contactus](https://www.thermofisher.com/contactus)

For product FAQs, protocols, training courses, and webinars, go to [thermofisher.com/technicalresources](https://www.thermofisher.com/technicalresources)



An example arrangement of five Applied Biosystems instruments receiving the MDTV service simultaneously.

Find out more at [thermofisher.com/thermalcyclerservice](https://www.thermofisher.com/thermalcyclerservice)

# High-performance PCR plastics for optimal PCR results

## MicroAmp PCR plastics

Applied Biosystems™ PCR plastics have been designed and validated to work with our thermal cyclers for more than 25 years. That's why they are engineer approved to enable optimal PCR performance.

### Applied Biosystems™ MicroAmp™ PCR plastics are:

- Validated on Applied Biosystems thermal cyclers for optimal fit and performance
- Designed for optimal heat transfer with thin-walled polypropylene wells
- Designed to reduce cross-contamination with raised well rims for effective sealing



### Applied Biosystems™ MicroAmp™ TriFlex 3 x 32-Well PCR Reaction Plate and adhesive film

- One plate that can be used as separated, individual 32-well plates or as a full 96-well plate
- Easy-to-tear dual side tabs
- Alphanumeric labeling on each plate segment
- Precut Applied Biosystems™ MicroAmp™ 32-Well Clear Adhesive Film

### Unique, high-performance features of Applied Biosystems™ MicroAmp™ EnduraPlate™ plastic consumables

#### Easy visual organization

5 color choices

#### Easy-to-read well identification text

Black text for excellent contrast

#### No warping, even after thermal cycling


Polycarbonate (hard shell) for mechanical stability and flatness

#### Snug fit to thermal blocks

Thin-wall polypropylene for excellent mechanical fit and heat transfer

#### Available in common formats

96-well standard, 96-well Fast, 384-well, 5-piece sample packs, 20 and 500 packs



Constructed to ANSI/SBS standard

Certified DNA-, RNase-, and PCR inhibitor-free

Compatible with and optimized for performance on Applied Biosystems instruments



### Options for every format and throughput need

Choose from tubes, tube strips, plates, adhesive film, and accessories for any throughput need. MicroAmp EnduraPlate plastics offer a solution for experiments that require special handling, such as automated or high-throughput workflows, and an even greater degree of durability for use with multi-instrument experiments.



MicroAmp 8-tube strip with attached domed caps; optical flat caps are also available for use with real-time PCR (qPCR).

The Applied Biosystems™ MicroAmp™ 8-Tube Strip with attached optical or domed caps offers a combination of features to help prevent cross-contamination, pipetting errors, and sample identification errors in your PCR and real-time PCR applications.

- Attached caps that open and close independently of each other
- Etched A–H letter labeling for individual tubes and caps
- Dual side tabs for strip labeling
- Graduated 20  $\mu$ L measuring markers on every tube
- Available exclusively for Applied Biosystems thermal cyclers to enable optimal PCR results


Find out more at [thermofisher.com/pcrplastics](https://thermofisher.com/pcrplastics)



### Did you know?

Proper plate sealing helps reduce evaporation and well-to-well contamination.

1. Remove the backing of the Applied Biosystems™ adhesive film.
2. Align the adhesive film so as to cover all wells while placing on the plate.
3. Rub the flat edge of the applicator along the long edge (length) of the plate, then along the short edge (width). Finally, rub the applicator between all the wells and around the outside edges of the plate using small back-and-forth motions to form a complete seal.



# Which PCR plastic fits your needs?

Find the PCR plastic format with the throughput and features for your application

Use for:	Small-scale experiments with a few samples	Routine experiments	Automation	Laboratory use
	Single tubes, strips, caps, adhesive film, and accessories	MicroAmp optical microplates	MicroAmp EnduraPlate optical microplates	MicroAmp EnduraPlate optical microplates GPLE*
<b>Formats</b>	<ul style="list-style-type: none"> <li>• Single tubes</li> <li>• Single tubes with caps</li> <li>• 8-strip tubes with caps</li> <li>• 12-strip caps</li> <li>• Adhesive film/plate seal</li> </ul>	<ul style="list-style-type: none"> <li>• 32-well</li> <li>• 48-well Fast</li> <li>• 96-well</li> <li>• 96-well Fast</li> <li>• 384-well</li> </ul>	<ul style="list-style-type: none"> <li>• 96-well</li> <li>• 96-well Fast</li> <li>• 384-well</li> </ul>	<ul style="list-style-type: none"> <li>• 96-well</li> <li>• 96-well Fast</li> <li>• 384-well</li> </ul>
<b>DNA-, RNase-, PCR inhibitor-free</b>	Yes	Yes	Yes	Yes
<b>ANSI/SBS standard dimension color</b>	Clear, or mixed packs containing red, orange, blue, and green	Clear	Single-color packs (red, blue, green, yellow, or clear) and 5-plate sampler (1 of each color)	Clear
<b>Instrument compatibility</b>	Use our plastics selection tool	Use our plastics selection tool	Use our plastics selection tool	Use our plastics selection tool
<b>Barcode available</b>	No	Yes (1 or 2 sides)	Yes (3 sides)	Yes (3 sides)
<b>Multiple applications</b>	No	No	Yes	Yes
<b>Optical compatibility</b>	Yes (applicable for optical version)	Yes	Yes	Yes
<b>Use</b>	Research use only	Research use only	Research use only	For laboratory use**

\*General purpose laboratory equipment (GPLE); for laboratory use.

\*\* Lot-based contamination test with Certificate of Analysis.




## Did you know?

Need high-quality PCR plastics for non-Applied Biosystems instruments?

Visit [thermofisher.com/thermoscientificplastics](https://thermofisher.com/thermoscientificplastics) for a wide range of Thermo Scientific™ PCR plastics.

Custom and OEM plastics for PCR and qPCR are available.

Learn more at [thermofisher.com/oemplastics](https://thermofisher.com/oemplastics)





## Quickly find the plastics and accessories you need for your instrument

Applied Biosystems MicroAmp Product	Cat. No.	3 x 32-well	96-well			96-well Fast	384-well		Genetic analyzers	
		ProFlex	ProFlex, SimpliAmp, VeritiPro, Veriti, MiniAmp Plus, MiniAmp	2720	9700	Veriti	ProFlex, VeritiPro	9700	310	3130, 3130xl, 3500, 3500xL, 3730, 3730xl
<b>96-well, 0.2 mL reaction plates</b>										
Optical 96-Well plate	N8010560, 4316813		•		• •					•
Optical 96-Well plate with barcode	4306737, 4326659		•		• •					•
Optical 96-Well plate with barcode and optical caps	403012		•		• •					
Optical 96-Well plate with barcode and optical adhesive films	4314320		•		• •					
EnduraPlate Optical 96-Well Clear plate	A36924		•		• •					•
EnduraPlate Optical 96-Well Clear plate with barcode*	4483354, 4483352		•		• •					•
TriFlex 3 x 32-Well PCR plate	A32810, A32811	•**	•		• •					
<b>96-well, 0.1 mL reaction plates</b>										
Fast Optical 96-Well plate	4346907					•			•	•
Fast Optical 96-Well plate with barcode	4346906, 4366932					•			•	•
EnduraPlate Optical 96-Well Fast Clear plate	A36930					•			•	•
EnduraPlate Optical 96-Well Fast Clear plate with barcode*	4483485, 4483494					•			•	•
<b>384-well reaction plates</b>										
Optical 384-Well plate	4343370						•	•		•
Optical 384-Well plate with barcode	4309849, 4326270, 4343814						•	•		•
EnduraPlate Optical 384-Well Clear plate	A36931						•	•		•
EnduraPlate Optical 384-Well Clear plate with barcode*	4483285, 4483273						•	•		•
<b>Strip tubes and caps</b>										
Fast 8-Tube Strip, 0.1 mL	4358293					•				
Optical 8-Tube Strip with Attached Optical Caps, 0.2 mL	A30588	•	•		• •					
8-Tube Strip with Attached Domed Caps, 0.2 mL	A30589	•	•		• •					
8-Tube Strip, 0.2 mL*	N8010580	•	•		• •				•	
Optical 8-Tube Strip, 0.2 mL	4316567	•	•		• •					
8-Cap Strip*, clear	N8010535, N8011535	•	•		• •	•				
Optical 8-Cap Strip	4323032	•	•		• •	•				
12-Cap Strip*	N8010534, N8011534	•	•		• •	•				
<b>Single tubes</b>										
Fast Reaction Tube with Cap, 0.1 mL	4358297, 4358293					•				
Reaction Tube with Cap, 0.2 mL*	N8010540, N8010612, N8011540	•	•		• •					
Reaction Tube without Cap, 0.2 mL*	N8010533, N8011533	•	•		• •					
Optical Tube without Cap, 0.2 mL	N8010933	•	•		• •					
<b>Seals and covers</b>										
Clear Adhesive Film	4306311		•		•	•	•	•		
Optical Adhesive Film	4360954, 4311971		•		•	•	•	•		
96-Well Full Plate Cover	N8010550				• •					
32-Well Clear Adhesive Film	A32812	•**	•		• •					
<b>Accessories</b>										
Splash-Free 96-Well Base	4312063		•		• •	•				
96-Well Support Base	4379590		•		• •	•				•
96-Well Base	N8010531		•		• •					

\* Multiple colors are available.

\*\* Do not use Applied Biosystems™ MicroAmp™ 3 x 32-Well Retainer (Cat. No. 4481669).

**Note:** Experiments using one or two 8-tube strips with attached caps require blank tube strips to balance lid pressure on the block or the use of the Applied Biosystems™ MicroAmp™ 96-Well Tray/Retainer Set (Cat. No. 4381850)—bottom part of tray only. For use with 96-well block of ProFlex, SimpliAmp, and VeritiPro thermal cyclers.

Visit our online plastics selection guide at [thermofisher.com/findplastics](https://thermofisher.com/findplastics)

# Superior cDNA synthesis for any application

## SuperScript IV Reverse Transcriptase

With more than 50,000 citations, reviews, and publications, Invitrogen™ SuperScript™ reverse transcriptases are among the most trusted and widely used products for cDNA synthesis. Invitrogen™ SuperScript™ IV Reverse Transcriptase is the latest enzyme in the portfolio, engineered to deliver superior cDNA synthesis performance with even the most challenging RNA samples.

- **Super-efficient**—up to 100x higher cDNA yields than with other reverse transcriptase enzymes
- **Super-sensitive**— $C_t$  values reduced by as many as eight cycles for RT-qPCR
- **Super-robust**—high thermostability and processivity for superior cDNA synthesis, even from degraded or inhibitor-containing RNA samples
- **Super-fast**—10 min cDNA synthesis protocol



### Select the right format for your application

We offer a comprehensive portfolio of enzymes and kits within the SuperScript IV family to suit your research needs. Start with the selection guide below to find the best format for common cDNA synthesis applications.

Would you like to have ...	... the ability to optimize reaction components and conditions?	... a complete kit with all cDNA synthesis reaction components?	... ultimate convenience and minimal pipetting steps for RT-PCR?	... the ability to skip the RNA extraction process and synthesize cDNA directly from mammalian cell lysate?
<b>Product format</b>	Stand-alone enzyme	First-strand cDNA synthesis kit	One-Step RT-PCR kit	Direct RT kit
<b>Recommended product</b>	SuperScript IV Reverse Transcriptase	SuperScript IV First-Strand Synthesis System	SuperScript IV One-Step RT-PCR System	SuperScript IV CellsDirect cDNA Synthesis Kit
<b>Applications</b>	RT-PCR, RT-qPCR, sequence detection, gene expression analysis, transcript variant detection, cloning, cDNA library construction, RACE, RNA-Seq	RT-PCR, RT-qPCR, sequence detection, gene expression analysis, transcript variant detection, cloning, cDNA library construction, RACE, RNA-Seq	RT-PCR, sequence detection, cloning, genotyping, high-throughput analysis	RT-PCR, RT-qPCR
<b>Input total RNA</b>	1 pg–5 μg	1 pg–5 μg	0.01 pg–1 μg	1–10,000 cells
<b>Optimal reaction temperature</b>	50–55°C	50–55°C	50–55°C	50–55°C
<b>Reaction time</b>	10 min	10 min	10 min	10 min
<b>cDNA synthesis with challenging or degraded RNA</b>	Yes	Yes	Yes	Yes

Order Invitrogen™ Oligo(dT) or random primers for cDNA synthesis at [thermofisher.com/rtpimers](https://thermofisher.com/rtpimers)

Find out more at [thermofisher.com/ssiv](https://thermofisher.com/ssiv)

# Complete kit with flexible priming options

## SuperScript IV First-Strand Synthesis System

The Invitrogen™ SuperScript™ IV First-Strand Synthesis System is optimized for synthesis of first-strand cDNA from purified poly(A)<sup>+</sup> or total RNA. The kit contains all components needed for reverse transcription, plus an additional control gene and primers, providing the flexibility to customize the reaction setup to fit the needs of your application.

Find out more at [thermofisher.com/ssiv-firststrand](http://thermofisher.com/ssiv-firststrand)



### Conducting a single-cell low-RNA input study?

Invitrogen™ SuperScript™ IV Single Cell/Low Input cDNA PreAmp Kit combines Invitrogen™ SuperScript™ IV and Invitrogen™ Platinum™ SuperFi™ enzymes along with template switching technology, which allows generation of high-quality full-length cDNA with high yields, sensitivity, and accuracy for RNA sequencing or qPCR. Find out more at [thermofisher.com/ssivpreamp](http://thermofisher.com/ssivpreamp)

# Enabling faster, more efficient RT-PCR

## SuperScript IV One-Step RT-PCR System

Even with challenging RNA samples, you can get more efficient results faster and easier than with any other RT-PCR reagent. The Invitrogen™ SuperScript™ IV One-Step RT-PCR System combines high-processivity SuperScript IV Reverse Transcriptase and high-fidelity Invitrogen™ Platinum™ SuperFi™ DNA Polymerase to provide superior one-step RT-PCR performance.

Find out more at [thermofisher.com/ssiv-onestep](http://thermofisher.com/ssiv-onestep)

- **Two-phase hot-start activation mechanism**—for high specificity, improved yields, and easy room temperature setup
- **Superior sensitivity and speed**—down to 0.01 pg of RNA, target length up to 13.8 kb, and the fastest one-step RT-PCR protocol
- **Reliable target detection**—even in RNA samples with suboptimal purity
- **Fast and easy gDNA removal**—for superior accuracy and confidence in your results

### Doing RT-qPCR?

Invitrogen™ SuperScript™ IV VILO™ Master Mix offers exceptional cDNA synthesis for RT-qPCR applications, while maintaining superior linearity across the broadest range of input RNA. Learn more when you flip the brochure to the real-time PCR side, or visit [thermofisher.com/4vilo](http://thermofisher.com/4vilo)

### Ordering information

Product	Quantity	Cat. No.
SuperScript IV Reverse Transcriptase	2,000 units	18090010
	10,000 units	18090050
	4 x 10,000 units	18090200
SuperScript IV First-Strand Synthesis System	50 reactions	18091050
	200 reactions	18091200
SuperScript IV First-Strand Synthesis System with ezDNase Enzyme	50 reactions	18091150
	200 reactions	18091300
SuperScript IV One-Step RT-PCR System	25 reactions	12594025
	100 reactions	12594100
SuperScript IV One-Step RT-PCR System with ezDNase Enzyme	50 reactions	12595025
	100 reactions	12595100
SuperScript IV VILO Master Mix	50 reactions	11756050
SuperScript IV CellsDirect cDNA Synthesis Kit	50 reactions	11750150
SuperScript IV Single Cell/Low Input cDNA PreAmp Kit	48 reactions	11752048

# DNA polymerases for PCR

A comprehensive portfolio of PCR enzymes and master mixes is available with the high performance and consistency you need. Start with the selection guide below to find the enzyme best suited for your PCR applications.

Visit our online selection guide at [thermofisher.com/pcrenzymes](http://thermofisher.com/pcrenzymes)



PCR type	High-fidelity PCR	Hot-start PCR	Hot-start PCR	Direct PCR
Do you need ...	... accurate DNA sequences?	... cleaner bands or to detect low-abundance targets?	... a chemical hot start?	... direct amplification without genomic DNA purification?
Recommended DNA polymerase	Platinum SuperFi II DNA Polymerase	Platinum II <i>Taq</i> Hot-Start DNA Polymerase	AmpliTaq Gold 360 DNA Polymerase	Platinum Direct PCR Universal Master Mix
<b>Applications</b>				
Cloning and subcloning	•			
Site-directed mutagenesis	•			
GC-rich templates	•	•	•	•
Sanger sequencing templates	•	•	•	•
High-throughput PCR	•	•		•
Long-range PCR (20–40 kb)	•			
Genotyping	•	•	•	•
Amplification of samples with suboptimal purity	•	•		•
Colony PCR	•	•	•	•
Multiplex PCR	•	•	•	•
Fast PCR	•	•		•
<b>Technical specifications</b>				
Fidelity compared to <i>Taq</i> DNA polymerase	>300x	1x	1x	1x
Target length	Up to 40 kb*	Up to 5 kb	Up to 5 kb	Up to 8 kb
Hot-start modification	Antibody-mediated	Antibody-mediated	Chemical modification	Antibody-mediated
Speed	15–30 sec/kb	15 sec/kb	1 min/kb	20 sec/kb
Universal primer annealing	Yes	Yes	NA	Yes
Inhibitor tolerance	Yes	Yes	NA	Yes
Blunt or 3'-A end	Blunt	3'-A	3'-A	3'-A
Compatibility with TaqMan probes	NA	Yes	Yes	NA
<b>Formats</b>				
Master mix	Colorless/green**	Colorless/green**	Colorless	Green**
Stand-alone enzyme	Colorless	Colorless/green†	Colorless	NA

\* Depends on complexity of DNA templates; may require additional optimization of reaction conditions and primer design.

\*\* Direct gel loading with green buffer options.

† Green buffer available as separate item for use with stand-alone enzyme.



## Did you know?

Assembled PCR reactions that include Invitrogen™ Platinum™ SuperFi™ II or Invitrogen™ Platinum™ II *Taq* Hot-Start DNA Polymerase are stable for 24 hours at room temperature, enabling high-throughput applications with automation or long waits between runs.

# Ultimate accuracy and robustness

## Platinum SuperFi II DNA Polymerase

Platinum SuperFi II DNA Polymerase is a hot-start, engineered, proofreading DNA polymerase for PCR applications requiring the highest sequence accuracy. Its fidelity is >300x that of *Taq* DNA polymerase and its buffer is specially formulated for a universal annealing temperature of 60°C.



- **Exceptional fidelity**—>300x more accurate than *Taq* DNA polymerase
- **Simplified workflow**—no need for primer melting temperature ( $T_m$ ) calculation; enables cocycling of multiple PCR targets
- **Increased PCR success**—robust amplification of GC-rich targets, DNA of suboptimal purity, and long sequences
- **Invitrogen™ Platinum™ hot-start technology**—enables superior specificity, sensitivity, and yields; allows for room temperature reaction setup and benchtop stability

Find out more at [thermofisher.com/platinumsuperfi](http://thermofisher.com/platinumsuperfi)

# Enhanced PCR specificity and yield with fast cycling

## Platinum II *Taq* Hot-Start DNA Polymerase

Platinum II *Taq* Hot-Start DNA Polymerase is an enzyme engineered for simplicity, specificity, and speed. Its universal primer annealing feature simplifies optimization and allows cocycling of different PCR assays.

- **Universal primer annealing at 60°C**—reduces tedious optimization steps and enables cocycling of different PCR assays

- **Engineered *Taq* DNA polymerase**—allows 4x faster synthesis and successful amplification even in the presence of inhibitors
- **Platinum hot-start technology**—enables superior specificity, sensitivity, and yields; allows for room temperature reaction setup and benchtop stability

Find out more at [thermofisher.com/platinumitaq](http://thermofisher.com/platinumitaq)

### Ordering information

Product	Quantity*	Cat. No.
Platinum SuperFi II DNA Polymerase	100 reactions	12361010
Platinum SuperFi II PCR Master Mix	100 reactions	12368010
Platinum SuperFi II Green PCR Master Mix	100 reactions	12369010
Platinum II <i>Taq</i> Hot-Start DNA Polymerase	100 reactions	14966001
Platinum II Hot-Start PCR Master Mix (2X)	50 reactions	14000012
Platinum II Hot-Start Green PCR Master Mix (2X)	50 reactions	14001012

\* Additional product sizes available.



# Everything you need for reliable PCR

## PCR primers and dNTPs

### Primers

We offer high-quality Invitrogen™ custom DNA oligos for successful PCR. Our oligos synthesis services feature over 25 years of customer service and popular primer design tools, as well as an informative, easy-to-use ordering portal to meet your needs. Invitrogen custom DNA oligos offer:

- Purity that you require: desalting, cartridge, HPLC, or PAGE
- Scales that suit your research needs: 25 nmol, 50 nmol, 200 nmol, 1 µmol, and 10 µmol
- Product format you prefer: tube or plate
- Affordable price
- Enhanced bulk ordering



Find out more at [thermofisher.com/primers](https://thermofisher.com/primers)

The free Invitrogen™ OligoPerfect™ Primer Designer enables you to:

- **Gain confidence**—Primer3-based design algorithm
- **Speed up**—design primers for PCR or capillary electrophoresis (CE), up to 50 genes at the same time
- **Work smarter**—recognizes .txt and .fasta file types
- **Order with ease**—seamlessly integrates with the Invitrogen™ ordering portal
- **Store your data**—ability to save your projects

Try OligoPerfect Primer Designer at

[thermofisher.com/oligoperfect-designer](https://thermofisher.com/oligoperfect-designer)

### Complex and large-scale oligos

Easily customize your oligos based on the unique scale, modification, and purification requirements for your research.

- Milligram to >15-gram delivered yields
- Unmodified and modified DNA or RNA up to 100 bases
- ISO 9001– and ISO 13485–controlled process
- Custom formulations and packaging solutions
- OEM collaboration and commercial supply support available
- Dedicated technical project team



Our experienced technical team can help you with your scale-up from step one to launch.

Learn more at [thermofisher.com/largescaleoligos](https://thermofisher.com/largescaleoligos)

### dNTPs

Our dNTPs have been extensively tested and verified for use in a wide variety of molecular biology applications, including highly sensitive techniques such as RT-qPCR and next-generation sequencing.

Learn more at [thermofisher.com/dntp](https://thermofisher.com/dntp)



## Did you know?

We offer additional manufacturing services, including OEM scale-up and support?

Leverage our expertise and infrastructure. We can manufacture products with your labels and packaging.

Explore OEM and commercialization options for oligo synthesis at [thermofisher.com/oem](https://thermofisher.com/oem).

# Simplify nucleic acid electrophoresis

## Nucleic acid separation and analysis

### E-Gel precast agarose gels

Using precast agarose gels can simplify the nucleic acid electrophoresis workflow. Invitrogen™ E-Gel™ precast gels are ready-to-use, prestained agarose gels packaged in a disposable cassette. There are no gels to pour, buffers to make, staining or destaining steps to perform, or gel boxes to assemble. Just load your samples and start the run.

- **Ready to use**—precast with agarose and DNA stain packaged inside a disposable cassette
- **Time-saving**—ideal for customers looking to streamline workflows
- **Versatile**—designed for rapid analysis of PCR products, plasmid preparations, and restriction digests to improve genotyping workflows



**1** Load      **2** Run      **3** Analyze

Sample analysis in three simple steps—load, run, and analyze.



Approximate time to complete an E-Gel electrophoresis system workflow and a traditional DNA electrophoresis workflow

### Ordering information

Product	Quantity	Cat. No.
E-Gel Agarose Gels with SYBR Safe DNA Gel Stain, 1%	10 gels	A42100
E-Gel Agarose Gels with SYBR Safe DNA Gel Stain, 2%	10 gels	A42135
E-Gel EX Agarose Gels, 2%	10 gels	G401002
E-Gel 48 Agarose Gels with SYBR Safe DNA Gel Stain, 2%	8 gels	G820802
E-Gel 96 Agarose Gels with SYBR Safe DNA Gel Stain, 2%	8 gels	G720802
E-Gel CloneWell II Agarose Gels with SYBR Safe DNA Gel Stain, 0.8%	10 gels	G661818

Find out more at [thermofisher.com/egel](https://www.thermofisher.com/egel)



### Did you know?

E-Gel precast gels are available in a variety of formats for routine and high-throughput applications, with different stains and agarose percentages (0.8%, 1.2%, 2%, and 4%). To find the right gel for your needs, see the selection guide at [thermofisher.com/egelselection](https://www.thermofisher.com/egelselection)

# Electrophoresis elevated to empower

## E-Gel Power Snap Plus Electrophoresis System

The Invitrogen™ E-Gel™ Power Snap Plus Electrophoresis System combines the convenience of quick real-time nucleic acid analysis with high-quality visualization in both high- and low-throughput Invitrogen™ E-Gel™ products. It offers:

- **Fast and safe workflow**—uses Invitrogen™ E-Gel™ technology, Invitrogen™ SYBR™ Safe DNA Gel Stain, and a blue-light transilluminator
- **Convenience**—gel-running and documentation on the same small device
- **Easy analysis**—available on both the device and on the Thermo Fisher™ Connect Platform with Invitrogen™ iBright™ image analysis software



### Specifications

Invitrogen™ E-Gel™ Power Snap Plus Electrophoresis Device	
Dimensions (L x W x H)	29.1 x 22.0 x 10.1 cm
Weight	2.15 kg
Touchscreen LCD display	3.5 in. TFT module with capacitive touch
LED light	Blue LED center wavelength: 465 nm; full width, half max: 20 nm
Viewing surface dimensions	15.4 x 11.2 cm
Amber filter dimensions	16.0 x 13.0 cm
LED life span	50,000 hours
LED specification	Array of 39 high-power LEDs with emission wavelength at 465 ± 10 nm


Invitrogen™ E-Gel™ Power Snap Plus Camera	
Dimensions (L x W x H)	23.3 x 21.0 x 22.8 cm
Weight	2 kg
Touchscreen LCD display	8 in. TFT module with capacitive touch
Camera type	Color, complementary metal-oxide semiconductor (CMOS)
Gel-image resolution	4,208 x 3,120 pixels (13 MP), 8 bits
Dynamic range	69 dB
Image output	TIF, JPG, PNG
Lens f-number	2.8
Internal memory	64 GB SD card

### Ordering information

Product	Cat. No.*	Cat. No.**
E-Gel Power Snap Plus Electrophoresis Device	G9101	G9110
E-Gel Power Snap Plus Camera	G9201	G9200
E-Gel Power Snap Plus Electrophoresis System	G9301	G9311
E-Gel Power Snap Plus Electrophoresis System with E-Gel 48 1% agarose gels	G9341	G9331
E-Gel Power Snap Plus Electrophoresis System with E-Gel 48 2% agarose gels	G9342	G9332
E-Gel Power Snap Plus Electrophoresis System with E-Gel 96 1% agarose gels	G9391	G9381
E-Gel Power Snap Plus Electrophoresis System with E-Gel 96 2% agarose gels	G9392	G9382

\* North America/Europe/Middle East/Africa

\*\* Asia Pacific/Japan/Latin America/Greater China



### E-Gel precast agarose gels

E-Gel precast agarose gels in 48- and 96-well with the option of 11- and 12-well gel analysis can be used with the Invitrogen E-Gel Power Snap Plus Electrophoresis System, and are offered with easy-to-use SYBR Safe stain and at percentages suitable for your lab's needs.

Find out more at [thermofisher.com/powersnapplus](http://thermofisher.com/powersnapplus)

# Fully integrated routine electrophoresis

## E-Gel Power Snap Electrophoresis System

The Invitrogen™ E-Gel™ Power Snap Electrophoresis System combines rapid, real-time, low-throughput nucleic acid analysis with high-resolution image capture for ultimate convenience. The integrated design helps reduce workflow time and accelerate discovery.

### Features of the E-Gel Power Snap Electrophoresis System include:

- **Faster analysis**—go from sample loading to image capture in as few as 15 minutes
- **Simple operation**—intuitive user interface with large touchscreen and integrated operating system
- **Safer workflow**—minimize handling of hazardous chemicals when used with E-Gel precast gel cassettes



### Specifications

Invitrogen™ E-Gel™ Power Snap Electrophoresis Device	
Dimensions (L x W x H)	24 x 13 x 7 cm
Weight	1 kg
Touchscreen LCD display	7.7 x 4.4 cm
LED light	Blue LED; center wavelength: 465 nm; full width, half max.: 20 nm
Viewing surface dimensions	9 x 11 cm
Amber filter dimensions	8.6 x 10.5 cm
LED life span	50,000 hours
LED specification	Array of 12 high-power LEDs with emission wavelength at 465 ± 10 nm

Invitrogen™ E-Gel™ Power Snap Camera	
Dimensions (L x W x H)	26 x 13 x 15 cm
Weight	1 kg
Touchscreen LCD display	11.5 x 8.6 cm
Camera type	Color, CMOS
Gel-image resolution	1,600 x 1,944 pixels (3 MP), 8 bits
Dynamic range	68 dB
Image output	TIF (grayscale) and JPG (color) file formats
Lens f-number	2.8
Internal memory	32 GB

### Ordering information

Product	Cat. No.
E-Gel Power Snap Electrophoresis Device	G8100
E-Gel Power Snap Electrophoresis Camera	G8200
E-Gel Power Snap Electrophoresis System	G8300
E-Gel Power Snap Electrophoresis System Starter Kit, E-Gel EX gels (1%)	G8341ST
E-Gel Power Snap Electrophoresis System Starter Kit, E-Gel EX gels (2%)	G8342ST
E-Gel Power Snap Electrophoresis Device Starter Kit, E-Gel SizeSelect II gels (2%)	G8162ST

 To learn more, go to [thermofisher.com/powersnap](https://www.thermofisher.com/powersnap)



### Did you know?

Samples can be run up to twice as fast with E-Gel precast gels compared to conventional handcast gels. E-Gel precast gels are available in a variety of formats: from low to high throughput, with SYBR Safe stain, and with agarose percentages suitable for either general-purpose or high-resolution separations.

# Pour your own agarose gels

## Stains and reagents for electrophoresis

### Pour-your-own electrophoresis reagents

Choosing the right tools for nucleic acid electrophoresis can significantly improve and accelerate results, enabling you to address downstream applications sooner.

### DNA stains

Detection of nucleic acid samples in gels can be improved using fluorescent dyes that are safer and/or more sensitive than ethidium bromide. The Invitrogen™ SYBR™ Safe, SYBR™ Green I, and SYBR™ Gold stains are less toxic and/or more sensitive with lower background fluorescence compared to conventional ethidium bromide stain.

Find out more at [thermofisher.com/stains](http://thermofisher.com/stains)

### UltraPure reagents for electrophoresis

Invitrogen™ UltraPure™ reagents are specifically formulated to meet your nucleic acid analysis and purification needs. UltraPure agarose and reagents are made from highly purified biochemicals for maximum reliability and superior performance.

Find out more at [thermofisher.com/ultrapure](http://thermofisher.com/ultrapure)



### Ordering information

Product	Quantity	Cat. No.
SYBR Safe DNA Gel Stain	400 µL	S33102
SYBR Gold Nucleic Acid Gel Stain	500 µL	S11494
UltraPure Agarose	100 g	16500100
TrackIt 1 Kb Plus DNA Ladder	100 applications	10488085
UltraPure TAE Buffer, 10X	4 L	15558026
UltraPure DNase/RNase-Free Distilled Water	10 x 500 mL	10977023
UltraPure Ethidium Bromide, 10 mg/mL	10 mL	15585011



# Specialty enzymes for molecular diagnostics development and commercial supply

Whether you need flexibility in assay design or custom functional testing, we can help accelerate your nucleic acid-based assay development with our innovative DNA-free and lyo-ready enzymes.

## DNA-free PCR enzymes\*

We've pioneered the use of single-use systems (SUS) to manufacture DNA-free PCR enzymes. The advantages include:

- Manufactured in a completely closed system using dedicated or single-use equipment
- Verified free of contaminating DNA from host, human operator, and environment
- Produced in an ISO 13485-certified facility for high standards of quality

\* Invitrogen™ Platinum™ Taq DNA Polymerase, DNA-free, is available as a catalog product.

Learn more at [thermofisher.com/dna-free](https://thermofisher.com/dna-free)

## Lyo-ready DNA polymerases and reverse transcriptases

From Platinum Taq DNA polymerases to SuperScript reverse transcriptases, we offer the largest selection of lyophilization-compatible enzymes, providing:

- The same functional enzyme performance as with the conventional format
- Tailor-made solutions for your specific applications, including custom formulation
- Higher confidence in results with low residual DNA contamination (human and bacterial)

Learn more at [thermofisher.com/lyo-ready](https://thermofisher.com/lyo-ready)

## Molecular biology education

Have you ever had questions about PCR but didn't want to ask? Find answers in our online education hub for molecular biology that features free technical content, fun videos, educational webinars, and application notes. Learn about PCR, reverse transcription, and more—at your own pace and on your own time—to amplify your research to the next level.

Find out more at [thermofisher.com/molbioschool](https://thermofisher.com/molbioschool)



[Tips and considerations](#)



[Technical videos](#)



[Educational webinars](#)



[Application notes](#)



## Mobile app: PCR Quest lab game

Test your PCR knowledge with our lab game, PCR Quest, where you travel from lab to lab crushing the world's toughest diseases. Download the game at [thermofisher.com/pcrquest](https://thermofisher.com/pcrquest)



Turn the brochure over to learn  
about our real-time PCR solutions.

Find out more at [thermofisher.com/pcrworkflow](https://thermofisher.com/pcrworkflow)

**For Research Use Only. Not for use in diagnostic procedures.** © 2019–2021 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. TaqMan is a trademark of Roche Molecular Systems, Inc., used under permission and license. AmpliTaq Gold is a trademark of Roche Molecular Systems, Inc., used under permission and license. Apple is a trademark of Apple Inc. Google is a trademark of Google Inc. **COL35052 1222**



# Real-time PCR solutions

More efficiency to take your research further


## Discover more efficient real-time PCR (qPCR)

The family of most-cited Applied Biosystems™ qPCR solutions delivers the caliber of results you need, faster and easier than you could have imagined. Free yourself from extra steps, preventable rework, and other inefficiencies so you can focus on making meaningful discoveries.

Experience the difference the efficient workflow of Applied Biosystems qPCR solutions can make to your research success:

- Connected instruments incorporate smart technologies to enable hands-free operation and remote support, while a suite of cloud-based analysis software provides exceptional flexibility.
- Predesigned Applied Biosystems™ TaqMan™ Assays provide everything you need, including Applied Biosystems™ TaqMan™ probe and PCR primer sets, formulated to work right out of the box. No additional design, optimization, or post-run analysis steps are required.
- Complement your assay with an Applied Biosystems™ TaqMan™ master mix and experience the gold standard in qPCR performance. We also offer Applied Biosystems™ SYBR™ dye chemistry with advanced features, like dye tracking, that help prevent pipetting errors.

**More choices. More flexibility. More efficiency.**  
**Experience more from your qPCR.**





# Contents

Which qPCR instrument fits your needs?	4
Real-time PCR applications	6
Real-time PCR software	7
QuantStudio 3 and 5 Real-Time PCR Systems	8
QuantStudio 6 Pro and 7 Pro Real-Time PCR Systems	10
QuantStudio 12K Flex Real-Time PCR System	12
Orbitor RS2 Microplate Mover	14
<hr/>	
TaqMan Array Plates with RFID	15
<hr/>	
MicroAmp qPCR plastics	16
<hr/>	
Which qPCR plastic fits your needs?	18
<hr/>	
TaqMan and SYBR Green chemistries for qPCR	20
<hr/>	
TaqMan Assays to accelerate your research	21
<hr/>	
TaqMan and SYBR Green master mixes	23
<hr/>	
SuperScript IV VILO Master Mix	25
<hr/>	
SuperScript IV CellsDirect cDNA Synthesis Kit	26
<hr/>	
Instrument Connect app	26
<hr/>	
Digital PCR (dPCR)	27
QuantStudio Absolute Q Digital PCR System	28
Instrument services and support	30
<hr/>	



# Which qPCR instrument fits your needs?




## QuantStudio real-time PCR systems

	Ultimate simplicity QuantStudio 3 system	Total control QuantStudio 5 system
		
Colors	4 colors	5 or 6 colors (21 filter combinations)
Available formats*	96-well (0.1 mL block) 96-well (0.2 mL block)	96-well (0.1 mL block) 96-well (0.2 mL block) 384-well
Dimensions (H x W x D)	40 x 27 x 50 cm	40 x 27 x 50 cm
Block change	Fixed	Fixed
VeriFlex temperature control	Yes, 3 zones	Yes, 6 zones (96-well blocks only)
Throughput	Medium	Medium
21 CFR Part 11–enabled	Security	Security, auditing, e-signature package
Touchscreen	Yes, interactive	Yes, interactive
Cloud-enabled	Yes	Yes
Automation compatible	NA	NA
Key applications	<ul style="list-style-type: none"> <li>• Gene expression</li> <li>• miRNA profiling</li> <li>• SNP genotyping</li> <li>• Copy number variation</li> <li>• Protein thermal shift</li> <li>• High-resolution melt</li> <li>• Pathogen detection</li> </ul>	<ul style="list-style-type: none"> <li>• Gene expression</li> <li>• miRNA profiling</li> <li>• SNP genotyping</li> <li>• Copy number variation</li> <li>• Protein thermal shift</li> <li>• High-resolution melt</li> <li>• Pathogen detection</li> </ul>

\* Some blocks may not be available.

Some instruments are also available in a diagnostic format. Learn more at [thermofisher.com/qsdX](https://thermofisher.com/qsdX)

 For information on other Applied Biosystems™ PCR instruments, please visit [thermofisher.com/qpcrinstruments](https://thermofisher.com/qpcrinstruments)

Smart and connected QuantStudio 6 Pro system*	Smarter productivity QuantStudio 7 Pro system*	Maximum productivity QuantStudio 12K Flex system
		
5 colors	6 colors (21 filter combinations)	6 colors (21 filter combinations)
96-well (0.1 mL block) 96-well (0.2 mL block) 384-well	96-well (0.1 mL block) 96-well (0.2 mL block) 384-well TaqMan Array Card (384-well microfluidic card)	96-well (0.1 mL block) 96-well (0.2 mL block) 384-well TaqMan Array Card (384-well microfluidic card) OpenArray plates (3,072 through-holes)
54.7 x 33.8 x 52.5 cm	54.7 x 33.8 x 52.5 cm	73.8 x 50.1 x 66.1 cm
Interchangeable; no tools required	Interchangeable; no tools required	Interchangeable; no tools required
Yes, 3 zones (96-well blocks only)	Yes, 6 zones (96-well blocks only)	NA
High	High	Very high
Upgradeable	Security, auditing, e-signature package	Optional security, auditing, e-signature package
Yes, interactive	Yes, interactive	Yes
Yes	Yes	No
Upgradeable	Yes	Yes
<ul style="list-style-type: none"> <li>• Gene expression</li> <li>• miRNA profiling</li> <li>• SNP genotyping</li> <li>• Copy number variation</li> <li>• Protein thermal shift</li> <li>• High-resolution melt</li> <li>• Pathogen detection</li> <li>• Pharmacogenomics</li> </ul>	<ul style="list-style-type: none"> <li>• Gene expression</li> <li>• miRNA profiling</li> <li>• SNP genotyping</li> <li>• Copy number variation</li> <li>• Protein thermal shift</li> <li>• High-resolution melt</li> <li>• Pathogen detection</li> <li>• Pharmacogenomics</li> </ul>	<ul style="list-style-type: none"> <li>• Gene expression</li> <li>• miRNA profiling</li> <li>• SNP genotyping</li> <li>• Copy number variation</li> <li>• Protein thermal shift</li> <li>• High-resolution melt</li> <li>• Pathogen detection</li> <li>• Pharmacogenomics</li> <li>• Growing menu of qualified solutions</li> </ul>

 = cloud-enabled instrument

## Real-time PCR applications

Real-time PCR is used for sensitive, specific detection and quantification of nucleic acid targets. We have developed powerful assay design algorithms, optimized master mixes, intuitive data analysis software, and flexible instrumentation to help harness the power of qPCR across a rich and diverse set of applications. Explore our robust solutions for your qPCR-based research.

### Infectious disease research

See our growing catalog of sensitive, specific qPCR probe and primer sets to study microbe biology and pathogenicity, host response, and disease pathology.

### Food pathogen detection

Detect multiple bacteria in the same run, including *Salmonella*, *Campylobacter*, *E. coli* O157:H7, and *Listeria monocytogenes*.

### Waterborne pathogen detection

Designed to detect and monitor waterborne pathogens in recreational and drinking water supplies.

### Pharmaceutical analytics

Designed to detect mycoplasmas, viruses, and residual host cell contamination for pharmaceutical, cosmetics, and personal care product manufacturing.

### Qualified solutions

A growing menu of new and valuable content for you to use on your high-throughput qPCR systems, including pharmacogenomics, vaginal microbiota, and *CFTR* mutation analysis.

### Stem cell research

Solutions for analyzing stem cells, determining stemness, and studying gene regulation and translation in stem cells.

### Pharmacogenomics research

Predesigned TaqMan Assays for more than 175 ADME and CYP targets, including >95% of ADME core markers and a warfarin metabolism panel.

### Oncology and genetic disease research

Enabling robust, reliable detection and quantitation of markers for cancer and genetic diseases.

### Plant sciences and agricultural biotechnology

Instruments, reagents, and kits designed for plant researchers that enable remarkable agricultural discoveries—from improved crops that feed more people to sustainable biofuels.

Other key applications include gene expression, genotyping, and sequencing.

# Real-time PCR software

A suite of analysis modules for the QuantStudio instrument family is available on the Thermo Fisher™ Connect Platform, enabling users to access and analyze their data anytime and anywhere. Just create and sign in to your account to use the Connect Platform for secure, cloud-based data storage, scientific analysis apps, and peer collaboration tools. Use the asset management apps and web tools to schedule time on your lab's instruments via your mobile device.

## [thermofisher.com/connect](https://thermofisher.com/connect)

Applied Biosystems qPCR analysis modules are a set of modules that provide an online toolkit for the analysis of qPCR data. The software takes advantage of the Connect Platform to provide highly versatile analysis tools that are flexible, fast, and easy to use and facilitate functional understanding of qPCR and related data. These modules allow users to combine over 100 genotyping, expression, or qPCR experiments into a project and analyze the data within minutes. Featured modules include:

- Design and Analysis
- High-Resolution Melt (HRM) Analysis
- Presence/Absence Analysis
- Standard Curve
- Genotyping
- hPSC Scorecard Analysis
- Relative Quantification



### Design and Analysis

The Design and Analysis app offers the ability to create, edit, and analyze qPCR instrument files.



### High-Resolution Melt (HRM) Analysis

This HRM Analysis app is designed for post-PCR analysis to identify variation in nucleic acid sequences. The method is based on detecting small differences in PCR melting (dissociation) curves. It is enabled by high-brightness, dsDNA-binding dyes used in conjunction with qPCR instrumentation that has precise temperature ramp control, advanced data capture capabilities, and access to software designed specifically for HRM analysis.



### Presence/Absence Analysis

The Presence/Absence Analysis app analyzes and interprets qPCR data or post-read data to determine if a specific target sequence is present in a sample or not. It provides an easy-to-view presence/absence result in a plate grid view.



### Standard Curve

The Standard Curve app offers reliable quantification of unknown quantities of genes and enables importing of standard curves from other experiments, providing analysis flexibility.



### Genotyping

The Genotyping app includes improved visuals and integrated traces of allelic discrimination plots to allow thorough QC of SNP assays to accurately reflect the true signals versus background noise.



### hPSC Scorecard Analysis

The hPSC Scorecard Analysis app scores gene expression profiles generated using the Applied Biosystems™ TaqMan™ hPSC Scorecard™ Panel compared to a reference set of well-characterized pluripotent stem cell (PSC) lines. It provides raw data as well as box plots, heat maps, and correlation plots for export.



### Relative Quantification

The Relative Quantification app allows fast and powerful gene expression analysis with enhanced visual capabilities for relative quantification, including integrated correlation and volcano and cluster analysis with the ability to drill down to amplification plots.

Learn more at [thermofisher.com/qpcranalysisapps](https://thermofisher.com/qpcranalysisapps)

## Intuitive and easy to use for all levels of experience

### QuantStudio 3 and 5 Real-Time PCR Systems

The Applied Biosystems™ QuantStudio™ 3 and QuantStudio™ 5 Real-Time PCR Systems provide intuitive touchscreen usability, allowing you to stay connected to your data easily. They're designed for both new and experienced users who need simple and affordable instruments without compromising performance or quality.

#### Get a premium instrument at an affordable price

**Access, analyze, and share data anytime and anywhere with the Connect Platform**—cloud-enabled services that support remote monitoring of your runs in real time, analyze sophisticated datasets in minutes, securely store data, and share results online with colleagues across institutions and around the world

**Obtain results you can trust**—detect differences in target quantity as small as 1.5-fold in singleplex reactions, and obtain 10 logarithmic units of linear dynamic range


**Multiplex with ease**—up to six excitation and six emission filters offer 21 different color combinations, allowing a broad range of detection chemistries and maximum multiplexing

**Helps save valuable time**—QuantStudio 3 and 5 systems offer three or six independent temperature zones for flexibility to run multiple experiments simultaneously; fast thermal cycling enables results in less than 30 minutes

**Get up and running quickly**—instruments are factory-calibrated for accuracy, quick installation, and immediate use; preoptimized protocol templates help minimize training for new users, and the included SmartStart Orientation provides basic qPCR training and setup for both the Connect Platform and the Instrument Management tool

**Maximize benchtop space**—the compact instrument can be configured stand-alone or with a computer



 = cloud-enabled instrument

Find out more at [thermofisher.com/quantstudio3-5](https://thermofisher.com/quantstudio3-5)



## Specifications

	QuantStudio 3 system	QuantStudio 5 system
Sample capacity (wells)	96	96 or 384
Reaction volume	0.1 mL block: 10–30 µL 0.2 mL block: 10–100 µL	96-well, 0.1 mL block: 10–30 µL 96-well, 0.2 mL block: 10–100 µL 384-well block: 5–20 µL
Footprint (H x W x D)	40 x 27 x 50 cm	
Excitation source	Bright white LED	
Optical detection	4 coupled filters	96-well: 6 decoupled filters 384-well: 5 coupled filters
Excitation/detection range	450–600 nm/500–640 nm	96-well: 450–680 nm/500–730 nm 384-well: 450–650 nm/500–700 nm
Multiplexing	Up to 4 targets	96-well: up to 6 targets 384-well: up to 5 targets
2D barcode reading	Optional	
Heating/cooling method	Peltier	
Temperature zone function	3 VeriFlex zones	96-well: 6 VeriFlex zones 384-well: NA
Max. block ramp rate	0.2 mL block: 6.5°C/sec 0.1 mL block: 9.0°C/sec	0.2 mL block: 6.5°C/sec 0.1 mL block: 9.0°C/sec 384-well block: 6.0°C/sec
Average sample ramp rate	3.66°C/sec	
Temperature uniformity	0.4°C	
Temperature accuracy	0.25°C	
Run time	<30 min runs	96-well block: <30 min runs 384-well block: <35 min runs
Dye compatibility (name)	FAM/SYBR Green, VIC/JOE/HEX/TET, ABY/NED/TAMRA/Cy3, JUN, ROX/Texas Red	FAM/SYBR Green, VIC/JOE/HEX/TET, ABY/NED/TAMRA/Cy3, JUN, ROX/Texas Red, Mustang Purple, Cy5/LIZ, Cy5.5
Chemistry capabilities	Fast/standard	
Features to assist with 21 CFR Part 11 compliance	No	Yes, with no additional fees
Detection sensitivity	1 copy	
Sensitivity	Detect differences as small as 1.5-fold in target quantities in singleplex reactions	

## Ordering information

Product	Cat. No.	Instrument + 1-year extended warranty with AB Assurance Cat. No.*
QuantStudio 3 Real-Time PCR System (96-well, 0.1 mL block)**	A28136	A33777
QuantStudio 3 Real-Time PCR System (96-well, 0.2 mL block)**	A28137	A33779
QuantStudio 5 Real-Time PCR System (96-well, 0.1 mL block)**	A28138	A33619
QuantStudio 5 Real-Time PCR System (96-well, 0.2 mL block)**	A28139	A33624
QuantStudio 5 Real-Time PCR System (384-well block)**	A28140	A33628

\* Includes SmartStart Orientation.

\*\* Does not include computer. Additional Cat. Nos. are available that include laptop or desktop computer.

## Recommended plastics

96-well block	Cat. No.	96-well block Fast	Cat. No.	384-well block	Cat. No.
MicroAmp EnduraPlate Optical 96-Well Clear Reaction Plate with Barcode	4483354	MicroAmp EnduraPlate Optical 96-Well Fast Clear Reaction Plate with Barcode	4483485	MicroAmp EnduraPlate Optical 384-Well Clear Reaction Plate with Barcode	4483285
MicroAmp Optical 96-Well Reaction Plate with Barcode	4306737	MicroAmp Fast Optical 96-Well Reaction Plate with Barcode, 0.1 mL	4346906	MicroAmp Optical 384-Well Reaction Plate with Barcode	4309849
MicroAmp Optical 8-Tube Strip with Attached Optical Caps, 0.2 mL	A30588	MicroAmp Fast 8-Tube Strip	4358293	MicroAmp Optical 384-Well Reaction Plate	4343370
MicroAmp Optical Adhesive Film	4360954	MicroAmp Optical 8-Cap Strip	4323032	MicroAmp Optical Adhesive Film	4360954
		MicroAmp Optical Adhesive Film	4360954		

# Smarter productivity and improved workflow

## QuantStudio 6 Pro and 7 Pro Real-Time PCR Systems

The Applied Biosystems™ QuantStudio™ 6 Pro and 7 Pro Real-Time PCR Systems are designed to save space on the lab bench, increase format flexibility, and support quick and simple block changes. These systems also deliver a smart workflow experience with advanced features that are designed for efficiency and productivity.

### Personalized

- Automatically load your settings and plate setup
- Simply log in with facial authentication, no passwords required
- SmartStart Orientation, which includes one day of on-site training that covers basic instrument operation and maintenance, and a choice of hands-on application training (including reagents) on either absolute quantification, gene expression, or genotyping

### Efficient

- Minimize hands-on time with voice commands facilitating hands-free operations
- Eliminate manual steps to obtain plate layout, protocol, and assay information on the instrument when using Applied Biosystems™ TaqMan™ Array Plates with RFID (see [page 15](#) for more information)
- Maximize uptime with the push-button Smart Help feature to access technical support scientists for faster, more efficient troubleshooting and resolution
- Reduce downtime using the Smart Remote Support collaboration tool with real-time video and audio, and desktop support to guide you through repairs remotely

### Simple

- Streamlined workflow directly from touchscreen
- Simple, tool-free block changes



\* Automation is available for the QuantStudio 7 Pro system and is available for the QuantStudio 6 Pro system after upgrade.

### Productive

- Access data anytime and anywhere with cloud-enabled services from the Connect Platform
- Built-in help videos for commonly used activities
- Increased ergonomics due to larger touchscreen, improved screen angle, and optional Bluetooth™ keyboard

The QuantStudio 7 Pro system also features automation compatibility and a security, auditing, and e-signature (SAE) package.

Find out more at [thermofisher.com/quantstudiopro](https://thermofisher.com/quantstudiopro)

## Specifications

	QuantStudio 6 Pro system	QuantStudio 7 Pro system
Sample capacity (wells)	96, 384	96, 384, TaqMan Array Card
Reaction volume	0.1 mL block: 10–30 µL 0.2 mL block: 10–100 µL 384-well block: 5–20 µL	0.1 mL block: 10–30 µL 0.2 mL block: 10–100 µL 384-well block: 5–20 µL TaqMan Array Card: ~1 µL
Footprint (H x W x D)	54,7 x 33,8 x 52,5 cm	
Excitation source	Bright white LED	
Filter or color combinations	5	21
Multiplexing	5 targets	6 targets
Excitation/emission range	450–680 nm/500–730 nm	
2D barcode reading	Via USB connection	
Heating/cooling method	Peltier	
Temperature zone function	3 independent temperature zones	6 independent temperature zones
Maximum ramp rate	6,5°C/sec	
Average sample ramp rate	3,66°C/sec	
Temperature uniformity	0,4°C	
Temperature range	4–99,9°C	
Temperature accuracy	0,25°C	
Run time	Less than 30 min	Less than 30 min
Dye compatibility (name)	FAM/SYBR Green, VIC/JOE/TET/HEX, TAMRA/NED/ABY, ROX/JUN/Texas Red, LIZ/Cy5/Mustang Purple, CY5.5/Alexa Fluor dye	
Features to assist with 21 CFR Part 11 compliance	Upgradeable	Yes
Detection sensitivity	1 copy	
Sensitivity	Detect differences as small as 1.5-fold in target quantities in singleplex reactions	

## Ordering information

Product	Cat. No.	Instrument + extended warranty package Cat. No.*
QuantStudio 6 Pro Real-Time PCR System, 96-well, 0.1 mL block	A43160	A47200
QuantStudio 6 Pro Real-Time PCR System, 96-well, 0.1 mL block, laptop	A43167	A47202
QuantStudio 6 Pro Real-Time PCR System, 96-well, 0.1 mL block, desktop	A43181	A47204
QuantStudio 6 Pro Real-Time PCR System, 96-well, 0.2 mL block	A43159	A44288
QuantStudio 6 Pro Real-Time PCR System, 96-well, 0.2 mL block, laptop	A43166	A44290
QuantStudio 6 Pro Real-Time PCR System, 96-well, 0.2 mL block, desktop	A43180	A44292
QuantStudio 6 Pro Real-Time PCR System, 384-well	A43161	A45582
QuantStudio 6 Pro Real-Time PCR System, 384-well, laptop	A43168	A45584
QuantStudio 6 Pro Real-Time PCR System, 384-well, desktop	A43182	A45586
QuantStudio 7 Pro Real-Time PCR System, 96-well, 0.1 mL block	A43163	A47201
QuantStudio 7 Pro Real-Time PCR System, 96-well, 0.1 mL block, laptop	A43170	A47203
QuantStudio 7 Pro Real-Time PCR System, 96-well, 0.1 mL block, desktop	A43184	A47205
QuantStudio 7 Pro Real-Time PCR System, 96-well, 0.2 mL block	A43162	A44289
QuantStudio 7 Pro Real-Time PCR System, 96-well, 0.2 mL block, laptop	A43169	A44291
QuantStudio 7 Pro Real-Time PCR System, 96-well, 0.2 mL block, desktop	A43183	A44293
QuantStudio 7 Pro Real-Time PCR System, 384-well	A43164	A4558383
QuantStudio 7 Pro Real-Time PCR System, 384-well, laptop	A43171	A45585585
QuantStudio 7 Pro Real-Time PCR System, 384-well, desktop	A43185	A4558745587
QuantStudio 7 Pro Real-Time PCR System, TaqMan Array Card block	A43165	A47651
QuantStudio 7 Pro Real-Time PCR System, TaqMan Array Card block, laptop	A43172	A47653
QuantStudio 7 Pro Real-Time PCR System, TaqMan Array Card block, desktop	A43186	A47655

\* The extended warranty package includes the instrument, SmartStart Orientation, and a 1-year AB Assurance Plan with 1 planned maintenance (PM) visit. Packages are not available in all countries. Contact your sales representative for more information.

Recommended consumables can be found on [page 15](#) (RFID-tagged plates), [page 16](#) (plastics), [page 21](#) (assays), [page 23](#) (master mixes), and [page 25](#) (reverse transcriptases).

# Maximum productivity with minimum effort

## QuantStudio 12K Flex Real-Time PCR System

This one instrument enables multiple users to conduct a wide range of experiments, from low- to high-throughput sample processing and for virtually any PCR application, such as:

- Drug discovery
- Pharmacogenomics research
- MicroRNA profiling
- Agriculture molecular testing
- *CFTR* mutation analysis
- Vaginal microbiota research

### Miniaturization at a lower cost

Applied Biosystems™ OpenArray™ technology is a broadly applicable nanoliter fluidics platform for low-volume solution-phase reactions, and enables lower reagent and assay costs, and rapid parallel processing.


### Superior throughput

The Applied Biosystems™ QuantStudio™ 12K Flex Real-Time PCR System can simultaneously run up to four 3,072-reaction Applied Biosystems™ OpenArray™ plates in about 4 hours. Up to 110,000 data points can be produced in an 8-hour day with the addition of a single Applied Biosystems™ ProFlex™ PCR System.

### Outstanding flexibility

Easily switch between five available thermal cycling blocks: OpenArray plate, TaqMan Array Card, 384-well, and standard or Fast 96-well blocks. Load the thermal cycling block in less than 1 minute—no tools required.

### OpenArray workflow



Find out more at [thermofisher.com/quantstudio12k](https://thermofisher.com/quantstudio12k)



= automation-ready instrument

### Increased data integrity and quality control—with results you can trust

The integrated sample tracking and SAE module assist you in supporting 21 CFR Part 11 compliance.

### Application spotlight

#### Pharmacogenomics research

Pharmacogenomics is the study of drug efficacy based on a subject's unique genomic composition. The QuantStudio 12K Flex system and OpenArray technology provide a simple, cost-effective, and fast workflow for the analysis of mutations and copy number variants associated with drug metabolism enzyme (DME) genes.

## Specifications

QuantStudio 12K Flex system		
Block change design	96-well (10–100 $\mu$ L reactions) Fast 96-well (15–30 $\mu$ L reactions) 384-well (5–20 $\mu$ L reactions) TaqMan Array Card (~1 $\mu$ L reactions) OpenArray plate (33 nL reactions)	
Block change design	Block change from front in less than 1 min; no tools or service call required	
Excitation source	Enhanced OptiFlex system, white LED	
Instrument control	Instrument touchscreen, networked computer, or attached computer	
Detection channels	<ul style="list-style-type: none"> <li>Decoupled: 6 emission, 6 excitation (96-well/Fast, 384-well, TaqMan blocks)</li> <li>Coupled: 4 emission, 4 excitation (OpenArray block)</li> </ul>	
21 CFR Part 11 enablement	Optional software module	
Dimensions (H x W x D)	73.8 x 50.1 x 66.1 cm	
Weight	70 kg	
Remote monitoring	Available to monitor up to 15 networked instruments simultaneously	
96-well, 96-well Fast, 384-well, TaqMan Array Card blocks		OpenArray block
Detection channels	Decoupled: 6 emission, 6 excitation	Coupled: 4 emission, 4 excitation
Well-to-well variability	$\pm 0.25^{\circ}\text{C}$	$\pm 0.75^{\circ}\text{C}$
Max. block ramp rate	$3.0^{\circ}\text{C}/\text{sec}$ (384-well)	$3.0^{\circ}\text{C}/\text{sec}$
Run time	30 min expected (Fast 96-well block) 35 min (384-well block, using Fast master mix)	2 hr (gene expression) 4 hr (genotyping)
Demonstrated sensitivity	To 1 copy	To 1 copy
Dynamic range	To 9 logarithmic units	To 7 logarithmic units
Resolution	As low as 1.5-fold change for singleplex reaction	As low as 2-fold change for singleplex reaction

## Ordering information

Product	Cat. No.	Instrument + 1-year extended warranty with AB Complete Cat. No.*
QuantStudio 12K Flex Real-Time PCR System, OpenArray block with AccuFill System, desktop configuration	4471090	4480621
QuantStudio 12K Flex Real-Time PCR System with OpenArray block without AccuFill System, desktop configuration	4472380	NA
QuantStudio 12K Flex Real-Time PCR System with TaqMan Array Card block, desktop configuration	4471089	4480622
QuantStudio 12K Flex Real-Time PCR System 384-well block, desktop configuration	4471134	4480623
QuantStudio 12K Flex Real-Time PCR System 96-well Fast block, desktop configuration	4471088	4480625
QuantStudio 12K Flex Real-Time PCR System 96-Well block, desktop configuration	4471087	4480631
96-Well Block Upgrade Kit	4453543	NA
96-Well Fast Block Upgrade Kit	4453544	NA
384-Well Block Upgrade Kit	4453545	NA
TaqMan Array Card Block Upgrade Kit	4453546	NA
OpenArray Block with AccuFill System	4471067	NA
QuantStudio 12K Flex AccuFill Upgrade Kit (For existing AccuFill system users)	4471022	NA

\* Includes SmartStart Orientation.

## Recommended plastics

96-well block	Cat. No.	96-well block Fast	Cat. No.	384-well block	Cat. No.
MicroAmp EnduraPlate Optical 96-Well Clear Reaction Plate with Barcode	4483354	MicroAmp EnduraPlate Optical 96-Well Fast Clear Reaction Plate with Barcode	4483485	MicroAmp EnduraPlate Optical 384-Well Clear Reaction Plate with Barcode	4483285
MicroAmp Optical 96-Well Reaction Plate with Barcode	4306737	MicroAmp Fast Optical 96-Well Reaction Plate with Barcode, 0.1 mL	4346906	MicroAmp Optical 384-Well Reaction Plate with Barcode	4309849
MicroAmp Optical 8-Tube Strip with Attached Optical Caps, 0.2 mL	A30588	MicroAmp Fast 8-Tube Strip	4358293	MicroAmp Optical 384-Well Reaction Plate	4343370
MicroAmp Optical Adhesive Film	4360954	MicroAmp Optical 8-Cap Strip	4323032	MicroAmp Optical Adhesive Film	4360954
		MicroAmp Optical Adhesive Film	4360954		

Multiple colors are available for most Cat. Nos.



# The perfect solution to automate a variety of workflows

## Orbitor RS2 Microplate Mover

The easy-to-use, flexible, and reliable Thermo Scientific™ Orbitor™ RS2 Microplate Mover is a collaborative benchtop mover that provides superior reliability and improved process efficiency. Its innovative bidirectional telescoping arm, coupled with its expansive 360° workspace, provides exceptional reach and precision, making the Orbitor RS2 Microplate Mover your trusted laboratory productivity partner.

### Designed with you in mind

- Integrated barcode reader enables sample tracking, barcode transmission, and inventory management
- Plate detection in the gripper helps to eliminate labware handling errors and reduces the risk of lost samples
- Variable force settings make it compatible with many plate types for a variety of applications



### Specifications

General specifications		Axis/joint specifications			
Power requirements	110–230 V, automatic switching	Axis	Range	Velocity	Acceleration
Operating environment	4–40°C, noncondensing	Z (height)	575 mm	750 mm/sec	2,250 mm/sec <sup>2</sup>
Safety compliances	CE, CSA	Radius (reach)	±406 mm	200 mm/sec	3,600 mm/sec <sup>2</sup>
Weight	25 kg	Theta (rotation)	Infinite	225 deg/sec	675 deg/sec <sup>2</sup>
Motion control	Closed-loop servo, brushless AC motor	Gripper stroke	51 mm		
Control	Directly from PC, requires Windows™ 7 or 10 operating system				
Storage options					
Random access hotel	8 plates per hotel, maximum capacity of 24 plates, best for 384-well plates				
Sequential access stack	40 plates per stack, maximum capacity of 80 plates, best for 384-well and 96-well (0.2 mL and 0.1 mL) plates				

### Ordering information

Product	Cat. No.
Orbitor RS2 Microplate Mover Bundle	A43780
Orbitor RS2 Microplate Storage, Sequential Access Stacks (3-pack)	A43781
Orbitor RS2 Microplate Storage, Random Access Hotels (3-pack)	A43782
Dual QuantStudio System–Orbitor RS2 Microplate Mover Integration Kit	A43779

Find out more at [thermofisher.com/orbitors2](https://thermofisher.com/orbitors2)

# One simple step gets you to your data faster


## TaqMan Array Plates with RFID

TaqMan Array Plates with RFID and QuantStudio 6 Pro and 7 Pro systems create an improved user experience with smart features that eliminate the manual steps currently necessary to set up the plate file for your qPCR instrument run. Just insert the RFID plate, and the QuantStudio 6 Pro or 7 Pro system captures the plate setup information, including:


- Plate layout—assay IDs in each well
- Plate type
- Expiration date
- Cat. No.
- Lot. No.
- Reaction volume per well
- Passive reference dye
- Thermal protocol



### Conventional workflow



### Smart qPCR workflow



### Ordering information

Description	No. of assays	No. of samples	Cat. No. (96-well 0.1 mL)	Cat. No. (96-well 0.2 mL)
TaqMan Array Plates with RFID	NA	NA	A43823	
TaqMan Array, 96-well Plate; format 8	8	12	4413263	4413266
TaqMan Array, 96-well Plate; format 16/16 Plus	16	6	4413261/4413262	4413264/4413265
TaqMan Array, 96-well Plate; format 32/32 Plus	32	3	4413259/4413260	4391528/4391529
TaqMan Array, 96-well Plate; format 48/48 Plus	48	2	4413257/4413258	4391526/4391527
TaqMan Array, 96-well Plate; format 96/96 Plus	96	1	4413255/4413256	4391524/4391525

Find out more at [thermofisher.com/taqman-rfid](https://thermofisher.com/taqman-rfid)

# High-performance qPCR plastics for optimal qPCR results

## MicroAmp qPCR plastics

Applied Biosystems™ PCR plastics have been designed and validated to work with our thermal cyclers for more than 25 years. That's why they are engineer approved to enable optimal PCR performance.

### Applied Biosystems™ MicroAmp™ qPCR plastics are:

- Validated on Applied Biosystems thermal cyclers for optimal fit and performance
- Designed to perform on all Applied Biosystems qPCR instruments
- Designed for optimal heat transfer with thin-walled polypropylene wells
- Designed to reduce cross-contamination with raised well rims for effective sealing



Plastics

## Unique, high-performance features of Applied Biosystems™ MicroAmp™ EnduraPlate™ plastic consumables

### Easy visual organization

5 color choices

### Easy-to-read well identification text

Black text for excellent contrast

### No warping, even after thermal cycling


Polycarbonate (hard shell) for mechanical stability and flatness

### Snug fit to thermal blocks

Thin-wall polypropylene for excellent mechanical fit and heat transfer

### Available in common formats

96-well standard, 96-well Fast, 384-well, 5-piece sample packs, 20 and 500 packs



Constructed to ANSI/SBS standard

Certified DNA-, RNase-, and PCR inhibitor-free

Compatible with and optimized for performance on Applied Biosystems instruments



### Options for every format and throughput need

Choose from tubes, tube strips, plates, adhesive film, and accessories for any throughput need. MicroAmp EnduraPlate plastics offer a solution for experiments that require special handling, such as automated or high-throughput workflows, and an even greater degree of durability for use with multi-instrument experiments.



The Applied Biosystems™ MicroAmp™ 8-Tube Strip with attached optical caps\* is optimally designed for precise real-time PCR with lid and tube labeling, dual end tabs, and 20  $\mu$ L graduation marks on each tube to prevent pipetting errors. The 8-tube strips fit in all 0.2 mL Applied Biosystems real-time PCR instruments.

\* MicroAmp 8-Tube Strip with attached domed caps is also available for PCR.

Find out more at [thermofisher.com/findplastics](https://thermofisher.com/findplastics)



### Did you know?

Proper plate sealing helps reduce evaporation and well-to-well contamination.

1. Remove the backing of the Applied Biosystems™ adhesive film.
2. Align the adhesive film so as to cover all wells while placing on the plate.
3. Rub the flat edge of the applicator along the long edge (length) of the plate, then along the short edge (width). Finally, rub the applicator between all the wells and around the outside edges of the plate using small back-and-forth motions to form a complete seal.



# Which qPCR plastic fits your needs?

Find the plastic format with the throughput and features for your application

	Small-scale experiments with a few samples	Routine experiments	Automation	Laboratory use
<b>Use for:</b>	Single tubes, strips, caps, adhesive film, and accessories	MicroAmp optical microplates	MicroAmp EnduraPlate optical microplates	MicroAmp EnduraPlate optical microplates GPLE*
<b>Formats</b>	<ul style="list-style-type: none"> <li>• Single tubes</li> <li>• Single tubes with caps</li> <li>• 8-strip tubes with caps</li> </ul>	<ul style="list-style-type: none"> <li>• 48-well Fast</li> <li>• 96-well</li> <li>• 96-well Fast</li> <li>• 384-well</li> </ul>	<ul style="list-style-type: none"> <li>• 96-well</li> <li>• 96-well Fast</li> <li>• 384-well</li> </ul>	<ul style="list-style-type: none"> <li>• 96-well</li> <li>• 96-well Fast</li> <li>• 384-well</li> </ul>
<b>DNA-, RNase-, PCR inhibitor-free</b>	Yes	Yes	Yes	Yes
<b>ANSI/SBS standard dimension color</b>	Clear	Clear	Single-color packs (red, blue, green, yellow, or clear) and 5-plate sampler (1 of each color)	Clear
<b>Instrument compatibility</b>	Use our plastics selection tool	Use our plastics selection tool	Use our plastics selection tool	Use our plastics selection tool
<b>Barcode available</b>	No	Yes (1 or 2 sides)	Yes (3 sides)	Yes (3 sides)
<b>Multiple applications</b>	No	No	Yes	Yes
<b>Optical compatibility</b>	Yes (applicable for optical version)	Yes	Yes	Yes
<b>Use</b>	Research use only	Research use only	Research use only	For laboratory use**

\* General purpose laboratory equipment (GPLE); for laboratory use.

\*\* Lot-based contamination test with Certificate of Analysis.




## Did you know?

Need high-quality PCR plastics for non-Applied Biosystems instruments?

Visit [thermofisher.com/thermoscientificplastics](https://thermofisher.com/thermoscientificplastics) for a wide range of Thermo Scientific™ PCR plastics.

Custom and OEM plastics for PCR and qPCR are available.

Learn more at [thermofisher.com/oemplastics](https://thermofisher.com/oemplastics)





## Quickly find the plastics and accessories you need for your instrument

Applied Biosystems MicroAmp Product	Cat. No.	48-well	96-well 0.2 mL			96-well 0.1 mL			384-well
		StepOne	7000	7300, 7500	QuantStudio 3, 5, 6 Flex, 7 Flex, 12K Flex, 6 Pro, 7 Pro; ViiA 7; 7900HT	StepOnePlus	7500	QuantStudio 3, 5, 6 Flex, 7 Flex, 12K Flex, 6 Pro, 7 Pro; ViiA 7; 7900HT	QuantStudio 5, 6 Flex, 7 Flex, 12K Flex, 6 Pro, 7 Pro; ViiA 7; 7900HT
<b>96-well, 0.2 mL reaction plates</b>									
Optical 96-Well plate	N8010560, 4316813		•	•	•				
Optical 96-Well plate with barcode	4306737, 4326659		•	•	•				
Optical 96-Well plate with barcode and optical caps	403012		•	•	•				
Optical 96-Well plate with barcode and optical adhesive films	4314320		•	•	•				
EnduraPlate Optical 96-Well Clear plate	A36924			•**	•				
EnduraPlate Optical 96-Well Clear plate with barcode*	4483354, 4483352			•**	•				
<b>96-well, 0.1 mL reaction plates</b>									
Fast Optical 96-Well plate	4346907					•	•	•	
Fast Optical 96-Well plate with Barcode	4346906, 4366932					•	•	•	
EnduraPlate Optical 96-Well Fast Clear plate	A36930					•	•	•	
EnduraPlate Optical 96-Well Fast Clear plate with barcode*	4483485, 4483494					•	•	•	
<b>384-well reaction plates</b>									
Optical 384-Well plate	4343370								•
Optical 384-Well plate with barcode	4309849, 4326270, 4343814								•
EnduraPlate Optical 384-Well Clear plate	A36931								•
EnduraPlate Optical 384-Well Clear plate with barcode*	4483285, 4483273								•
<b>48-well reaction plates</b>									
Fast Optical 48-Well Plate	4375816	•							
<b>Strip tubes and caps</b>									
Fast 8-Tube Strip, 0.1 mL	4358293	•				•	•	•	
Optical 8-Tube Strip with Attached Optical Caps, 0.2 mL	A30588		•	•	•				
Optical 8-Tube Strip, 0.2 mL	4316567		•	•	•				
Optical 8-Cap Strip	4323032	•	•	•	•	•	•	•	
<b>Single tubes and caps</b>									
Fast Reaction Tube with Cap, 0.1 mL	4358297	•				•		•	
Optical Tube without Cap, 0.2 mL	N8010933		•	•					
<b>Seals and covers</b>									
Optical Adhesive Film	4360954, 4311971		•	•	•	•	•	•	•
48-Well Optical Adhesive Film	4375323	•							
<b>Reaction trays</b>									
96-Well Tray/Retainer Set	403081		•						
Fast 48-Well Tray	4375282	•							
96-Well Tray for VeriFlex Block	4379983					•			
<b>Accessories</b>									
Splash-Free 96-Well Base	4312063		•	•	•	•	•	•	
96-Well Support Base	4379590		•	•	•	•	•	•	

\* Multiple colors are available.

\*\* Requires use of proper adapter, Cat. No. A24820.

**Note:** Experiments using one or two 8-tube strips with attached caps require blank tube strips to balance lid pressure on the block or the use of the Applied Biosystems™ MicroAmp™ 96-Well Tray/Retainer Set (Cat. No. 4381850)—bottom part of tray **only**. For use with 96-well block of Applied Biosystems™ 7000, 7300, 7500, and ViiA 7 systems, and QuantStudio 3/5/6/7/12K instruments.

Visit our online plastics selection guide at [thermofisher.com/findplastics](https://thermofisher.com/findplastics)

# TaqMan and SYBR Green chemistries for real-time PCR

## We offer two types of chemistries to detect PCR products using real-time PCR instruments:

- Applied Biosystems™ TaqMan™ chemistry (also known as fluorogenic 5' nuclease chemistry)
- Applied Biosystems™ SYBR™ Green I dye chemistry

With more than 21 million predesigned assays, including 2.8 million gene expression assays for more than 32 species, 5 assay formats, and >200,000 publications, TaqMan Assays represent the most trusted and comprehensive collection of qPCR assays available.

Find master mixes, reagents, and kits to power your applications.

	TaqMan chemistry—based detection	SYBR Green dye—based detection
<b>Chemistry overview</b>	Uses a fluorogenic probe to enable detection of a specific PCR product as it accumulates during PCR cycles	Uses SYBR Green I, or similar dye that binds to double-stranded DNA to detect the PCR product as it accumulates during PCR

	TaqMan Assays and reagents	SYBR Green reagents
<b>Specificity</b>	High	Low
<b>Sensitivity—low copy number</b>	High	Variable*
<b>Reproducibility</b>	High	Variable*
<b>Multiplexing</b>	Yes	No
<b>Predesigned assays</b>	Yes	No
<b>Custom assays</b>	Yes	No
<b>User design and optimization</b>	No	Yes
<b>Cost</b>	High	Low*
<b>Gene expression quantitation</b>	High	Low
<b>DNA quantitation</b>	Yes	Yes (pathogen detection)
<b>ChIP</b>	Yes	Yes
<b>SNP genotyping</b>	Yes	No
<b>MicroRNA</b>	Yes	No
<b>Copy number</b>	Yes	No
<b>Somatic mutation detection</b>	Yes	No
<b>Pathway analysis</b>	Yes	No

\* Depends on template quality, and primer design and optimization.

Chemistry	Application	Starting material	Recommended master mix
<b>5' nuclease (TaqMan Assays)</b>	Gene expression, DNA quantification	cDNA, gDNA	TaqMan Fast Advanced Master Mix
	Gene expression, RNA virus quantification	RNA	TaqMan Fast Virus 1-Step Master Mix
	Genotyping	gDNA	TaqPath ProAmp Master Mix
<b>SYBR Green dye</b>	Gene expression	cDNA	PowerTrack SYBR Green Master Mix

Find your assay at [thermofisher.com/taqman](https://thermofisher.com/taqman) or your master mix at [thermofisher.com/qpcrmastermixes](https://thermofisher.com/qpcrmastermixes)

# TaqMan Assays to accelerate your research

## Comprehensive, high-quality solutions for genetic analysis

TaqMan Assays are the industry-leading choice for 5' nuclease qPCR assays. They are cited in more publications than any other qPCR assay product and are considered the gold standard for quantitative genomic analysis. Backed by a performance guarantee,\* TaqMan Assays are consistently chosen as a proven solution to reliably provide fast and accurate results.

- **Specificity**—advanced primer/probe sequence selection criteria plus minor groove binder (MGB) probe enhancement deliver the specificity and reproducibility you need for confidence that your results are generated from amplification of the intended target and not from nonspecific dye binding or amplification of closely related genes or pseudogenes
- **Sensitivity**—the nonfluorescent quencher (NFQ) on Applied Biosystems™ TaqMan™ probes minimizes background, and intelligent PCR primer and probe design maximizes amplification efficiency; get better sensitivity and accuracy to reliably detect targets present at 10 or fewer copies
- **Reproducibility**—accurately reproduce results from well to well, day to day, and lab to lab, even across manufacturing lots
- **Proven technology**—referenced by more than 200,000 publications to date

**Guaranteed to perform for all your research needs\***



**TaqMan Assays cover a broad range of research for a wide variety of applications:**

Expression	Genetic variation
<ul style="list-style-type: none"> <li>• Gene expression</li> <li>• MicroRNA (miRNA) expression</li> <li>• Long noncoding RNA (lncRNA) expression</li> <li>• Fusion transcript detection</li> <li>• Protein expression</li> </ul>	<ul style="list-style-type: none"> <li>• Single-nucleotide polymorphism (SNP) genotyping</li> <li>• Drug metabolism enzyme genotyping</li> <li>• Copy number variation (CNV)</li> <li>• Rare somatic mutation detection</li> </ul>

\* Terms and conditions apply. See full details of the guarantee at [thermofisher.com/taqmanwarranty](https://thermofisher.com/taqmanwarranty)

### Applied Biosystems™ TaqMan™ Gene Expression Assays

- Designed to detect virtually any gene product, with more than 2.8 million pre-designed assays
- Best-coverage assays available to detect the highest number of transcript variants possible
- Available for 32 species and some microbial pathogens
- Multiple format options—single tubes, 96-well and 384-well plates, TaqMan Array Cards, and OpenArray plates
- Can't find a pre-designed assay to meet your needs? Use our Custom Assay Design Tool to design an assay based on your specific target sequence

### Applied Biosystems™ TaqMan™ Advanced miRNA Assays

- Detect targets with as few as 60 copies of input microRNA (miRNA) in the cDNA synthesis reaction
- One universal reverse transcription (RT) step for all miRNA assays
- cDNA can be archived for future miRNA studies
- Identify only mature miRNA and distinguish related highly homologous miRNAs with gold-standard TaqMan probe specificity
- Detect and quantify mature miRNA from as little as 1 pg of total RNA or 2 µL of purified plasma or serum

### Applied Biosystems™ TaqMan™ Noncoding RNA Assays

- For reliable detection and quantitation of noncoding transcripts longer than 200 nucleotides
- Designed to detect only noncoding transcript targets

### Applied Biosystems™ TaqMan™ Fusion Assays

- Quantify the expression level of gene fusions using qPCR
- Orthogonal validation method for confirming next-generation sequencing (NGS) results

### Applied Biosystems™ TaqMan™ SNP Genotyping Assays

- Comprise the world's largest predesigned collection, with more than 17 million assays
- Rigorous design pipeline with >90 parameters for optimal primer–probe combinations
- Functional quality testing on at least 20 gDNA samples for each assay
- Flexible formats to accommodate any number of targets and samples

### Applied Biosystems™ TaqMan™ Drug Metabolism Genotyping Assays

- Detect polymorphisms in high-value drug metabolism enzyme (DME) pharmacogenetics markers
- Cover specific SNP alleles, multinucleotide polymorphisms (MNPs), and insertions and deletions (indels)
- Targets derived from public databases, consortia, and published articles

### Applied Biosystems™ TaqMan™ Copy Number Assays

- Evaluate copy number of genomic DNA targets
- Easy to interpret—Applied Biosystems™ CopyCaller™ Software provides the calculated copy number and predicted copy number, along with confidence value and z-score quality metrics
- Fast and simple—setup to primary analysis in 3–4 hours

### Applied Biosystems™ TaqMan™ Mutation Detection Assays

- Powered by competitive allele-specific Applied Biosystems™ TaqMan™ PCR (castPCR™) technology
- Detect and measure somatic mutations in genes associated with cancer research
- Detect rare amounts of mutated DNA in a sample that contains large amounts of normal, wild-type DNA
- Compatible with different sample types, such as cell lines; formalin-fixed, paraffin-embedded (FFPE) tissue; and fresh-frozen tissue samples

### Custom TaqMan Probes

Design your own TaqMan primers and probes by choosing from dual-labeled Applied Biosystems™ TaqMan™ MGB, QSY™, or TAMRA™ custom probes. TaqMan MGB probes include an MGB moiety at the 3' end that increases the melting temperature ( $T_m$ ) of the probe and stabilizes probe–target hybrids. TaqMan MGB probes can be significantly shorter than traditional probes, providing better sequence discrimination and flexibility to accommodate more targets. TaqMan QSY probes are available for multiplexing three or more targets and are a seamless replacement for BHQ™ probes.

## Flexible formats

A variety of formats to meet your system and throughput needs

### TaqMan Assay formats


#### Configurations to fit your research goals

Are you analyzing hundreds (or thousands) of samples, and expression from a handful of genes? Or does your research involve a few samples that need to be analyzed for a long list of mRNA targets? No matter your experimental design, there's a TaqMan Assay format\* that fits your research needs.



#### Single tubes

- Low entry price
- Flexible for small projects
- Run on any qPCR instrument



#### 96- or 384-well plates

- Optimal for medium to large projects
- Balance flexibility with streamlined reaction setup
- Run on 96- or 384-well qPCR instrument



#### 384-well TaqMan Array Cards

- Low cost per reaction
- Optimal for medium to large projects
- Run on QuantStudio 7 Pro, 7 Flex, and 12K Flex systems\*\*

\* Also available on OpenArray plates for the QuantStudio 12K Flex system.

\*\* Also available on Applied Biosystems™ ViiA™ 7 and 7900HT systems.

# TaqMan and SYBR Green master mixes

TaqMan and Applied Biosystems™ SYBR™ Green master mixes have been designed to give you the reliable results you need in run times as short as 40 minutes. Master mixes are optimized for specific applications and contain everything needed for reliable qPCR—buffer, dNTPs, passive reference dye, thermostable hot-start DNA polymerase, and other components. Simply add your sample and assay, and start your reaction.



## Applied Biosystems™ TaqMan™ Fast Advanced Master Mix

Accurate, dependable gene expression quantification in less time

- Wide linear dynamic range even on fast cycling instruments (<40 min run time)
- 72-hour pre-PCR benchtop stability enables flexibility in automated workflows
- Compatible with single and multiplex assays

## Applied Biosystems™ TaqMan™ Fast Virus 1-Step Master Mix

Sensitive 1-step quantification of RNA viruses and transcripts, even with challenging samples

- Designed for reliable, high-sensitivity 1-step RT-qPCR to enhance virus and transcript quantification
- Formulated to handle common RT-PCR inhibitors found in blood, stool, and other difficult samples
- Single-tube 4X master mix format enables use of more sample for enhanced sensitivity
- Works in singleplex and multiplex, and with exogenous or endogenous internal controls

## Applied Biosystems™ TaqPath™ ProAmp™ and ProAmp™ Multiplex Master Mixes\*

Accurate SNP and CNV genotyping calls in both purified and crude samples

- Reproducible genotyping results even in the presence of inhibitors
- Two formulations that enable single- to four-target detection per reaction
- Compatible with buccal swab and blood crude lysates prepared with the Applied Biosystems™ DNA Extract All Reagents Kit
- 72-hour pre-PCR benchtop stability enables flexibility in automated workflows

## Applied Biosystems™ TaqPath™ BactoPure™ Microbial Detection Master Mix\*

Less background DNA for lower-level microbial detection

- Low-level detection of bacterial, fungal, mammalian, and viral DNA, as well as of antibiotic resistance markers
- Tolerance of PCR inhibitors typically found in both purified samples and crude lysates in biopharmaceutical, molecular diagnostic, and research applications
- Wide dynamic range\*\* (up to 8 orders of magnitude) enabling accurate detection from both low- and high-concentration samples

Explore Applied Biosystems™ QualTrak™ qPCR solutions for biopharma discovery at [thermofisher.com/qpcr/biopharma](https://thermofisher.com/qpcr/biopharma)

\* TaqPath ProAmp and TaqPath BactoPure master mixes are General Purpose Reagents for Laboratory Use.

\*\* Dynamic range is a property of both the assay and template concentration in the sample, as well as the formulation of the master mix; thus, individual results may vary.

Learn more about master mixes and how to use our selection tool at [thermofisher.com/qpcrmastermixes](https://thermofisher.com/qpcrmastermixes)



## Applied Biosystems™ PowerTrack™ SYBR™ Green Master Mix

### SYBR Green dye–based gene expression quantification with unrivaled specificity

- Helps reduce pipetting errors with a two-color tracking dye system
- Formulated for superior specificity and tight reproducibility of  $C_t$  values over a broad dynamic range
- Broad primer  $T_m$  and primer concentration compatibility allows flexibility in qPCR reaction setup with minimal optimization
- Compatible with Invitrogen™ SuperScript™ IV VILO™ Master Mix for fast, reproducible reverse transcription results

## Applied Biosystems™ PowerUp™ SYBR™ Green Master Mix

### Exceptional performance for even the most challenging qPCR applications

- High specificity with proprietary *Taq* DNA polymerase with dual-lock, hot-start mechanism
- Room-temperature stability of preassembled qPCR reactions for up to 72 hours
- Compatibility with standard or fast cycling for results in less than 50 minutes

### Ordering information (more sizes available)

Product	Quantity	Cat. No.
TaqMan Fast Advanced Master Mix	1 x 1 mL	4444556
	1 x 5 mL	4444557
	1 x 50 mL	4444558
	1 x 1 mL (no UNG)	A44359
	1 x 5 mL (no UNG)	A44360
TaqMan Fast Virus 1-Step Master Mix	1 x 1 mL	4444432
	5 x 1 mL	4444434
	1 x 10 mL	4444436
TaqPath ProAmp Master Mix	1 x 1 mL	A30865
	1 x 10 mL	A30866
	1 x 50 mL	A30867
TaqPath ProAmp Multiplex Master Mix	1 x 1 mL	A30868
	1 x 10 mL	A30869
TaqPath BactoPure Microbial Detection Master Mix	1 x 50 mL	A30870
	1 x 1 mL	A52699
	1 x 5 mL	A52700
PowerTrack SYBR Green Master Mix	1 x 50 mL	A52702
	1 mL	A46012
	5 mL	A46109
PowerUp SYBR Green Master Mix	50 mL	A46113
	1 mL	A25741
	5 mL	A25742
	50 mL	A25743

# Superior cDNA synthesis performance in RT-qPCR applications

## SuperScript IV VILO Master Mix

SuperScript IV VILO Master Mix is a first-strand cDNA synthesis reaction mix for two-step RT-qPCR. The master mix format elevates the trusted VILO technology (variable input, linear output) to the next level by combining further optimized buffer conditions with the highly processive and thermostable Invitrogen™ SuperScript™ IV Reverse Transcriptase. The master mix offers exceptional performance features while maintaining superior linearity across the broadest range of input RNA.

- **Super-efficient**— $C_t$  values earlier by an average of 2 cycles compared to other reverse transcription reagents, in a 10 min reaction
- **Super-strong**—reliable results even with degraded or inhibitor-containing RNA samples
- **Super-reliable**—improved RT-qPCR data reproducibility due to single-tube master mix format
- **Super-safe**—integrated, easy, and RNA-friendly genomic DNA removal

Find out more at [thermofisher.com/4vilo](https://thermofisher.com/4vilo)



**Doing other cDNA synthesis applications?**

The SuperScript IV Reverse Transcriptase portfolio of products is engineered to offer superior cDNA synthesis performance with even the most challenging RNA samples. The portfolio includes a stand-alone enzyme, first-strand cDNA synthesis kit, and one-step RT-PCR kit.

[thermofisher.com/ssiv](https://thermofisher.com/ssiv)

Reagents

### Ordering information

Product	Quantity	Cat. No.
SuperScript IV VILO Master Mix	50 reactions	11756050
	500 reactions	11756500
SuperScript IV VILO Master Mix with ezDNase Enzyme	50 reactions	11766050
	500 reactions	11766500



### Did you know?

SuperScript IV VILO Master Mix is available in a format with the novel dsDNA-specific Invitrogen™ ezDNase™ Enzyme, which offers convenient and fast genomic DNA removal from RNA samples to help ensure high confidence in RT-qPCR results.

# Why isolate RNA when you don't have to?

## SuperScript IV CellsDirect cDNA Synthesis Kit

The Invitrogen™ SuperScript™ IV CellsDirect™ cDNA Synthesis Kit is designed to synthesize first-strand cDNA directly from mammalian cell lysates without first isolating the RNA. With lysis and reverse transcription (RT) performed in the same tube, the resulting first-strand cDNA is ready to use in many downstream applications such as PCR and qPCR.



### Key benefits of this easy-to-use format

- **Time-saving**—save up to 2 hours in overall workflow time
- **Superior sensitivity**—easy detection of low-abundance targets
- **Direct convenience**—no need to isolate RNA
- **qPCR compatibility**—superior performance with PowerTrack SYBR Green Master Mix

Find out more at [thermofisher.com/cellsdirect](http://thermofisher.com/cellsdirect)

### Ordering information

Product	Cat. No.
<b>SuperScript IV CellsDirect cDNA Synthesis Kit</b>	
50 reactions	11750150
500 reactions	11750350
<b>SuperScript IV CellsDirect Lysis Reagents</b>	
500 reactions	11750550

## Connect with your instrument and achieve lab-life balance

### Instrument Connect app

Thermo Fisher Scientific has pioneered a way to connect you to your instruments, giving you real-time updates on your run and access to data as they come up. We connect you to your cloud-enabled instruments and benchtop devices anytime and anywhere using the Instrument Connect mobile app.

The Instrument Connect remote monitoring app allows you to stay connected to any of our cloud-enabled instruments, including the QuantStudio 3, 5, 6 Pro, and 7 Pro qPCR instruments, as well as endpoint PCR devices including Applied Biosystems™ ProFlex™, SimpliAmp™, and MiniAmp™ thermal cyclers.


With the Instrument Connect app, you can:

- Check the availability of your cloud- and network-connected device
- Monitor run progress
- View amplification plots in real time (available for QuantStudio 3, 5, 6 Pro, and 7 Pro qPCR instruments)
- View plots and filter by sample or target in real time
- Schedule an instrument

It's easy to get started.

Just download the Instrument Connect app from the Apple™ App Store or Google Play™ Store and log in to your account on the Connect Platform. You can view your connected instruments, monitor remaining time in your run, and view your amplification plots in real time.

Learn more about the Instrument Connect app at [thermofisher.com/connect](http://thermofisher.com/connect)





## Discover powerfully simple digital PCR

Digital PCR (dPCR) is quickly becoming a reliable, complementary technology to qPCR for precise nucleic acid quantification because it can deliver a high level of precision, sensitivity, and robust performance. dPCR overcomes common limitations of qPCR, such as the need for routine standard curves, low precision when measuring rare targets, and lack of sensitivity in high-background conditions.

Experience the difference that the efficient workflow of Applied Biosystems dPCR solutions can make to your research success:

- **Single instrument, single plate, simple workflow**—eliminate cumbersome reaction preparation steps from your dPCR workflow. The Applied Biosystems™ QuantStudio™ Absolute Q™ Digital PCR System requires only one hands-on step that takes less than five minutes to complete with minimal technical skill.
- **Fast time to answer**—avoid wasting precious time generating microchambers, moving plates, and reading samples. In 90 minutes, reagent digitization, thermal cycling, and data collection are integrated into a single system with no manual transfer steps required.
- **Industry-leading choice**—the Applied Biosystems brand is relied on for genetic analysis, offering robust assays powered by our unique design algorithms and years of experience. Our easy-to-use predesigned assay technology makes dPCR simple with Applied Biosystems™ Absolute Q™ digital PCR assays.

# dPCR that works how you need it, when you need it

## QuantStudio Absolute Q Digital PCR System



The QuantStudio Absolute Q Digital PCR System is a plate-based dPCR platform powered by proprietary microfluidic array plate (MAP) technology that enables all the necessary steps for dPCR—compartmentalizing, thermal cycling, and data acquisition—to be conducted on a single instrument. The dPCR workflow is identical to the qPCR workflow you are familiar with, improving ease of use, minimizing hands-on steps, and maximizing consistency.

The higher precision, sensitivity, and absolute nature of dPCR is ideal for research in:

- Rare-target detection, such as somatic mutation detection in oncology research
- Precise quantification of viral targets
- Pathogen detection and load determination
- Generation of references and standards
- Copy number variation

### Simple workflow

The simplified dPCR workflow requires one instrument, one consumable, and one hands-on step that takes less than 5 minutes to complete with minimal technical skill.

### Fast turnaround time

In 90 minutes, you can generate highly precise dPCR data in a fraction of the time required for typical dPCR workflows.

### Industry-leading consistency

MAP technology enables exceptional consistency in total microchambers analyzed per reaction.

### Maximum reagent efficiency

More than 95% of the input sample is analyzed per reaction compared to the 25–60% range from other dPCR platforms.

### Plate flexibility

Each 16-sample dPCR plate may be run one column at a time to preserve unused wells for future runs, minimizing waste.

### Absolute quantification

The >20,000 microchambers per sample generate data expressed as copies/ $\mu$ L, enabling high precision and sensitivity without the need for a standard curve.

### Application spotlight


#### Rare cancer mutation quantification in liquid biopsy

Liquid biopsies use cell-free DNA (cfDNA) that is derived from both normal and cancerous cells that have undergone apoptosis or necrosis and have released their DNA contents into their environment. The small fraction of cfDNA molecules originating from the tumor are called circulating tumor DNA (ctDNA), and these show cancer mutations unique to the tumor. dPCR is an established method for the fast, reliable, and accurate quantification of these cancer mutations that occur in low abundance. dPCR offers outstanding specificity and precision, making it a cost-effective and rapid research method for routine monitoring of ctDNA from liquid biopsy samples.


Preparation




Digitization



Amplification



Quantitation



Find out more at [thermofisher.com/absoluteq](https://thermofisher.com/absoluteq)



## Specifications

QuantStudio Absolute Q Digital PCR System	
Time to results	<2 hr
Colors	4 (endpoint detection)
Illumination	Rax, blue, phosphor green high-power LED
Detection channels	FAM, VIC, and HEX
Reader size (H x W x D)	54.0 x 60.0 x 62.0 cm
Weight	60 kg
QuantStudio Absolute Q Digital PCR MAP	
Samples per plate	16
Targets per plate	4
Microchambers per sample	20,480
Loading volume	9 µL
Sample waste/dead volume	<5%
Performance	
Dynamic range	5 logarithmic units
Precision at 95% confidence interval	±10%
Compatible chemistries	Absolute Q Digital PCR Assays and TaqMan Assays

## Ordering information

Product	Cat. No.
QuantStudio Absolute Q Digital PCR System, desktop, extended warranty, SmartStart Orientation	A53267
QuantStudio Absolute Q Digital PCR System, desktop	A52864
QuantStudio Absolute Q MAP16 Plate Kit and Master Mix	A53301
QuantStudio Absolute Q MAP16 Plate Kit	A52865
Absolute Q DNA Digital PCR Master Mix (5X)	A52490
QuantStudio Absolute Q Isolation Buffer	A52730

# Keep working wonders. We've got your back.

Count on superior services and support to keep your real-time and digital PCR instruments up and running



## Service plans help protect your investment

All new instruments are protected by a standard factory warranty. This warranty covers all costs for travel, labor, and parts for repairs. Extended coverage service plans are also available at the time of instrument purchase. Services at additional charge include:

- Service plans with guaranteed response times\*
- Priority technical support—priority phone access to our trained and certified technical and instrument support specialists
- 24-hour response time—expedited access to remote service engineers within 24 hours

Check your instrument's coverage status at [thermofisher.com/manage](https://www.thermofisher.com/manage)

## Digital service innovations

Our instrument service plans include digital innovations that help keep your instruments and your lab running smoothly. With pioneering on-demand tools and capabilities such as remote support using mobile augmented-reality (AR) technology, instrument-driven support, and on-demand instrument training, we're constantly looking ahead so your lab doesn't fall behind.

Explore our digital service innovations at [thermofisher.com/innovations](https://www.thermofisher.com/innovations)

## Instrument hardware qualification (IQ/OQ/PQ)

Our manufacturer-trained and certified field service engineers will conduct and document comprehensive tests, including software and hardware compatibility matrixing, component verification, and site requirements, to verify performance and provide reliable, audit-style documentation to meet your regulatory requirements. Installation qualification (IQ), operational qualification (OQ), and performance qualification (PQ) or instrument performance verification (IPV) services are recommended at installation, after moving the instrument, after software or hardware upgrades, and after planned maintenance or critical repairs.

Contact an instrument qualifications specialist at [thermofisher.com/iqoqpq](https://www.thermofisher.com/iqoqpq)

## Analytical validation consulting services

Why handle your analytical validation (AV) in-house, when using our AV consulting services can speed up the launch process by up to 75% and cut costs by up to half? We'll consult with you to develop and optimize your validation workflow, while providing data analysis support and template documentation to fully maximize your instrument and reagent investment.

Get a complete list of assays and panels we support, or contact an AV consulting services specialist at [thermofisher.com/av](https://www.thermofisher.com/av)

\* Guaranteed response times are dependent on type of service plan and geographic region.

## Service plans at a glance

	On-site service plans			Repair Center service plan*
	AB Platinum	AB Assurance**	AB Maintenance Plus	AB Repair Center Support Plus
Response time	Next business day <sup>†</sup>	2 business days <sup>†</sup>	3 business days <sup>†</sup>	
Planned maintenance	•	•	•	
Access to technical support (Mon–Fri, standard business hours)	24/7/365	•	•	•
Parts, labor, and travel	•	•	10% discount	•
Digital remote support	•	•	•	•
Qualification service	•	Available as add-on	Available as add-on	
Field application scientist (FAS) consultation	•	Available as add-on	Available as add-on	
Loaner instrument issued during repair (Repair Center plan only)				Available as add-on

\* Repair Center service plan available for QuantStudio 3 and 5 qPCR systems only.

\*\* AB Assurance service plan offered for QuantStudio Absolute Q Digital PCR System.

† Availability limited in some geographic areas.

‡ After receipt of purchase order.

### Education services

It can be difficult to prepare yourself for what's next while you're focused on the work you have now. Our professional, interactive training courses make it easier.

We offer a combination of virtual and in-person classroom instruction, and hands-on learning in your lab to match your schedule, budget, and learning preferences. Whichever course style you choose, you'll learn from one of our 300 highly skilled application scientists who are available to lead sessions online, at your location, or at one of our 12 training centers located worldwide.

Explore courses at [thermofisher.com/education/services](https://www.thermofisher.com/education/services)

### Technical support

If you have questions about product selection or use, assay or experimental design, data analysis, or troubleshooting, contact our team of technical support scientists or access our online product and application support tools.

### How to reach us

To find your local support or technical support team, go to [thermofisher.com/contactus](https://www.thermofisher.com/contactus)

For product FAQs, protocols, training courses, and webinars, go to [thermofisher.com/technicalresources](https://www.thermofisher.com/technicalresources)

Explore our services and support solutions at [thermofisher.com/instrumentservices](https://www.thermofisher.com/instrumentservices)



Turn the brochure over to learn  
about our PCR solutions

 Find out more at [thermofisher.com/qpcr](https://thermofisher.com/qpcr)

**For Research Use Only. Not for use in diagnostic procedures.** © 2019–2022 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. TaqMan is a trademark of Roche Molecular Systems, Inc., used under permission and license. AmpliTaq Gold is a trademark of Roche Molecular Systems, Inc., used under permission and license. Cy is a registered trademark of GE Healthcare. Apple is a trademark of Apple Inc. Google Play is a trademark of Google Inc. Bluetooth is a trademark of Bluetooth SIG, Inc. Windows is a trademark of Microsoft Corp. BHQ is a trademark of Biosearch Technologies Inc. **COL35052 1222**